


Extended antigen specific cellular responses by single-cell immune profiling

Abstract:
Evaluating the overall antigen-specific cellular immunity of an individual is important for understanding anti-tumor immunity and for development of personalized vaccines and immunotherapy. The dCODE Dextramer® (RIO) and dCODE Klickmer®, technology together with the BD Rhapsody™ Single-Cell Analysis workflow combine single-cell genomic profiling with antigen-specific recognition allowing deep phenotypic and genotypic analysis of antigen-specific T and B cells. However, to unveil the heterogeneity of the cellular immune responses seen in cancer and infectious disease many specificities need to be investigated to get the complete picture. Here we report how multiple samples can be assayed for antigen specific: T cells, MAIT, iNKT cells, and B cells responses, simultaneously in a single workflow.




Experimental setup:
We sample tagged 3 different cell samples individually, and composed a single sample in which we would demonstrate detection of antigen specific cells in the 5 different immune cell types: CD8+ cytotoxic T cells, CD4+ T helper cells, iNKT cells, MAIT cells, and B cells (Fig. 12).

The 3 samples and sample tag (SMK):
SMK1: A healthy pre-COVID-19 HPBMC sample with several known virus specific cytotoxic CD8 T cells.
SMK2: A CD4 enriched cell sample, stimulated and expanded in the presence of a EBV specific peptide.
SMK12: A normal HPBMC sample from a previously COVID-19 infected and vaccinated healthy Donor.
The cells were mixed in a ratio of app. 10(SMK1) to 1(SMK2) to 15(SMK12) respectively, then stained with a panel of dCODE Dextramer® (RIO) reagents (Table 1), and thereafter the Immune Discovery Panel (IDP), containing 30 BD® AbSeq assay (Table 2). dCODE Dextramer® positive cells (PE+) were sorted, and the cells were subjected to the BD Rhapsody™ System for full length TCR/BCR, VDJ sequencing, BD® AbSeq assay and dCODE Dextramer® (RIO) Library Preparation, and DNA sequencing. After sequencing of the libraries, the seq. data were processed through the Seven Bridges, BD Rhapsody™ Sequence Analysis Pipeline and subsequent analyzed using BD SeqGeq™ software package.

Data analysis
The data from the composed sample were analyzed in whole for its canonical markers and gated for detection of each of the antigen specific cell types (Fig. 2 a-e). A “canonical marker” parameter were generated containing CD3, CD4, CD8, CD19, and CD161, plotted against CD14, to remove monocytes from the downstream analysis. The CD14+ cells were gated for CD3+ and CD19+ to separate the T cells from the B cells, and the CD4+ and CD8+ positive T cells were subsequently gated from the CD3+ positive population. The CD14+ cells were furthermore gated for CD161+ and MR1(5-OP-RU)+ antigen specific MAIT cells, and for CD1d(α-GalSer)+ iNKT cells. Each of the CD4+, CD8+ and CD19+ populations were analyzed for the presence of antigen specific cells defined in the Dextramer® panel (Table 1, Fig. 2, f-g plots). The antigen specific cells were immunophenotyped using the IDP panel (Table 2) and defined as being, Naive, Naive like, Effector cells, Central memory cells and Effector memory cells (Fig. 4). The found antigen specific cells were furthermore defined by its full length clonotype (Fig. 6).

Fig 6: Results of Antigen-specific Cell Clonal Analysis:
Populations of antigen-specific CD4+ helper, and CD8+ cytotoxic T cells, Spike specific B cells, as well as iNKT and MAIT cells, were identified, and the full length VDJ sequences were analyzed. VDJ composition of the paired T cell and B cell receptors are shown.



171 cells were sequenced with paired TCR. Consistency with MAIT cells being “invariable”, we found low variability of the TCRs identified and the major part being TRAV V1-2*01, with low variability of the other domains, as been reported elsewhere.
*MR1-Restricted T Cells Are Unprecedented Cancer Fighters
Alessandro Vacchini, Andrew Chancellor, Julian Spagnuolo, Lucia Mori and Genaro De Libero*
Experimental Immunology, Department of Biomedicine, University of Basel and University Hospital Basel, Basel, Switzerland

6 cells were sequenced with paired TCR. Consistency with iNKT cells being “invariable” T cells, we found low variability of the TCRs identified. The VDJ alpha region represent only two one V and one J domain and a single α CDR3 seq. the β-VJ region, represent only one V domain, but 5 different J domains, with varying CDR3 sequences. Although low number of cells, it show the invariability of at least the α region, and the VDJ usage of both α and β regions.

9 cells of a single T cell clone were detected and sequenced with paired TCR. Consistency with EBV being an endogenous virus, that over time will converge to few antigen specific T cell clones.

The MHC Dextramer® panel generated 37 antigen specific CD8+ T cells. The sequencing of these cells was not optimal, but we could recognize the CDR3 clonotypes of the highest frequent T cell clone

We identified 15 COVID Spike specific B cells in this sample. The variable domains are significant different between the B-cell clones.

Fig 2: Gating main cell types, and identification of antigen specific cells




Fig 2: a-e gating of the 5 immune cell types in the composed sample; a) Exclusion of CD14+ monocytes, b) CD3+ T cells & CD19+ B cells, c) CD4+ and CD8+ T cells, d) CD161+MR1+MAIT cells, and e) CD3+CD1d+iNKT cells.
Lower f-h: Identification of the antigen specific T and B cells, overlaid (Blue) of the deconvoluted sample tag (single sample) antigen specific cells. Identifying the same cells in the composed sample, as in the sample tag deconvoluted samples (blue).

Fig 3: Use of internal neg. cell control to define antigen specific CD8 T cells




Fig.3: Sample tag (SMK) can be used to introduce an internal cell control for defining antigen specific responses. In this experiment we wanted to identify MHCII responses in an allele matched sample tagged donor (SMK1). We used a donor with no overlapping alleles, sample tagged with SMK12 as internal cellular control. CD8+ T cells were identified in each donor as a population, and the “gene expression / parameter density” of each was plotted in gene-view, CD8+SMK1 vs CD8+SMK12 expression, and then parameters were selected for either:
a) The 5 negative control dCODE MHCII Dextramer® reagents or b) The analytic dCODE MHCII Dextramer® virus specificities.
a) No differences are observed between CD8+ population of SMK1 and SMK12 for the negative control reagents.
b) 5 easily identifiable MHCII responses are “overexpressed” in SMK1 vs SMK12. Five MHCII specificities were found negative in SMK1 vs SMK12 plot. The same responses were found / not found when evaluating each MHCII specificity against CD8+ individually (data not shown).

Fig 4: Phenotyping antigen specific B-cells, CD8+ T cells and CD4+ T cells




Fig 4: Phenotyping antigen specific cells; B-Cells, CD8 T Cells, and CD4 T cells using BD® AbSeq assay antibodies. The antigen specific cells overlaid on the antigen negative cells, are shown.
a) Most of the identified COVID-19, Spike Specific B cells are IgM and IgD negative, indicating B cells type shifting toward mature blast cells. b) The Virus specific CD8 T cells differ from the main T cell population by being primarily CD45 and CCR7 negative, and CD27+, being Central memory cells. c) The EBV CD4 T cells (very few) are CD45, and CD27- having an effector Memory phenotype.

Fig 5: Phenotyping and placing the antigen specific cells on the cellular atlas




Fig 5: PCA was performed on all parameters, and the Principal components with a variance >1% were used to generate the t-SNE plot.
a) t-SNE with overlay of the major cell types defined by BD® AbSeq assay antibodies.
b) t-SNE with overlaid antigen specific responses, showing:
The dCODE MHCII Dextramer®+ cells are colocalized within the CD8+CD56- population, the COVID-19 SPIKE+ B-Cells within the CD19+B-Cell pop. The EBV class II positive cells are within the CD4+ populations. CD1d(α-GalSer)+ iNKT cells are mainly in the CD8+CD56+ population of NK-cells, and some within the CD4+ T cell population. MR1(5-OP-RU)+ cells, are all within the CD8+CD56+ T cell population and seems to be the driver in segregating the CD8+ cells into two populations of the t-SNE plots.

Summary:

We here demonstrate analysis of antigen-specific immune responses in 5 antigen-specific immune cell types: CD8+ cytotoxic T cells, CD4+ T helper cells, iNKT, and MAIT cells, as well as B cells, using a single workflow on a composed sample, in which each cellular component is sample tagged, individually.

Multiple virus specific responses were identified at the single cell level and cells characterized based on receptor recognition and full length VDJ clonotype sequence. Full length VDJ TCR/BCR parried clonotypes were identified for each of the responses identified.

The antigen responses are found with equal resolution and sensitivity, whether on the background of the whole composed sample or the SMK deconvoluted individual samples, indicating that mixing multiple sample is not decreasing the overall sensitivity of the assay.

Using sample tagging to introduce an internal negative control sample in the assay, enables identification of antigen specific T cells with high specificity, based on a “differential expression” analysis between the two sample tagged samples, showing that sample tagging can be an efficient tool where an internal cellular control can be introduced in the assay.

Furthermore, we enable phenotypic characterization of the antigen specific cells, using a 30-member immune response panel of BD® AbSeq antibodies (IDP), as well as placed the found antigen specific cells within each of the overall immune cell types, showing antigen specific cells were identified only in the cellular compartment that makes biological sense.

Table 1: dCODE Dextramer® (RIO) reagent panel.

Cell No.	dCODE Dextramer® (RIO) pool:	Antigen	Cell type
FA10002	HLA-DRB1*01:01/PSYRHRMTRALMQR (ELIP)	NC	CD4 T helper cell
FA10003	DRB1*01:01/TSYLNLRKGTALA (EBV)	EBV	CD4 T helper cell
WA5380	HLA-A*01:01/SLGGGLGY	NC	Cytotoxic T cell
WR466	HLA-A*02:01/LAPVHVAV	NC	Cytotoxic T cell
WF3231	HLA-A*24:02/AYSSAGASL	NC	Cytotoxic T cell
W13397	HLA-B*07:02/GPASAAGAL	NC	Cytotoxic T cell
WB2161	HLA-B*08:01/AAKRGAGAL	Flu	Cytotoxic T cell
WB2162	HLA-B*08:01/GELPVFTL	Flu	Cytotoxic T cell
WB2130	HLA-B*02:01/LGLTAVML	EBV	Cytotoxic T cell
WB2132	HLA-B*02:01/RLVAVML	EBV	Cytotoxic T cell
WB2133	HLA-B*02:01/VLETSVML	ChV	Cytotoxic T cell
WK2146	HLA-B*35:01/EPLPGQLTAY	EBV	Cytotoxic T cell
WK2145	HLA-B*35:01/HPVGEADYFYL	EBV	Cytotoxic T cell
WK2138	HLA-B*35:01/PSNVVHHY	ChV	Cytotoxic T cell
WA2131	HLA-A*01:01/VEHTDLY	Flu	Cytotoxic T cell
WA3410	HLA-A*01:01/CTELKLSY	ChV	MAIT cell
XB0802	HCID16 (α-GalCer)	Metabolit	iNKT cell
XB0801	HCID16 (α-GalCer)	NC	iNKT cell

Table 2: BD® AbSeq Immune Discovery Panel

Specificity	Clone	Oligo ID	Cell type / phenotype
CD3	IGCH1	AHS0231	canonical T cell marker
CD4	SA2	AHS0032	canonical T helper cells
CD8	SK1	AHS0228	canonical Cytotoxic T cells
CD13c	Bv4c	AHS0056	differentiation marker
CD14	MPI199	AHS0037	canonical Monocyte marker
CD16	308	AHS0053	Natural killer cell marker
CD19	3252C1	AHS0030	differentiation marker
CD25	2A3	AHS0026	T reg. and T cell activation marker
CD27	HW1231	AHS0025	differentiation marker
CD28	L293	AHS0138	T cell activation marker
CD38A	HL108	AHS0009	Naive cell marker
CD56	NCAM16	AHS0019	canonical Natural killer cells
CD62L	DR6G-56	AHS0049	differentiation marker (naive/central memory cells)
CD127	HL176-HEL	AHS0020	differentiation marker (Effector and memory cells)
CD134	ACT15	AHS0013	T cell activation marker
CD137	484	AHS0014	T cell activation marker
CD161	HP-3G10	AHS0005	MAIT cell marker & Natural killer cell marker
CD183 (CXCR3)	16C/CXCR3	AHS0031	T cell activation marker
CD183 (CXCR3)	318B8	AHS0039	T follicular helper (TFH) cells
CD186 (CXCR6)	13B 165	AHS0148	HIV coreceptor
CD186 (CXCR6)	118R	AHS0034	B cell activation marker
CD197 (CCR7)	Z-11-A	AHS0073	Naive differentiation marker
CD272	J168-540	AHS0052	Naive T cell marker
CD278	D929	AHS0012	T cell activation marker
CD279 (PD-1)	EH12-1	AHS0014	T cell exhaustion marker
CD325 (GITR)	V27-580	AHS0010	T cell activation marker
CD366 (TIM-3)	7D3	AHS0016	T cell exhaustion marker
HLA-DR	G46-G	AHS0035	T cells Differentiation / activation marker
IGD	L46-2	AHS0098	B cell Differentiation marker T cells
IGM	G20-127	AHS0198	B cell Differentiation marker T cells

This work was sponsored by BDbiosciences®