

MHC II dCODE Dextramer technology allows characterization of antigen-specificity, TCR clonotype and gene expression of single CD4+ T-cells

Bjarke Endel Hansen*, Sara Bursomanno, Kivin Jakobsen, Amir Ameri, Liselotte Brix


Immudex, Fruebjergvej 3, 2100 Denmark.
*Correspondance: beh@immudex.com

Background


T-cell mediated adaptive immunity to pathogens, tumor- and self-antigens in autoimmunity is dependent on the specific recognition of the peptide-MHC (pMHC) complex on the target cell by a unique T-Cell receptor (TCR). Identifying the specific T-cells, their phenotypic gene expression profiles and clonotypic TCR sequence is essential to understand the complexity of an immune response and to manipulate it for therapeutic benefits. We have developed an MHC multimer technology, dCODE™ Dextramer®, comprising a dextran backbone displaying multiple pMHCs and a unique DNA barcode coding for the displayed pMHC specificity. Detection of cells labelled with dCODE Dextramer is based on sequencing of the attached DNA barcode. The dCODE Dextramer is compatible with 10x Genomics Chromium platform for single cell analysis.

In this study, a panel of three different MHC II dCODE Dextramer reagents was used to identify antigen-specific CD4+ T-cells in a cell sample, together with their unique TCR clonotypes and gene expression profiles.

Two workflows were compared. In one, cells stained with dCODE Dextramer were directly used to perform single cell analysis, while in the other, stained cells were previously enriched for antigen specific T-cells by cell sorting.



Workflow and Results



A) The three DNA libraries were sequenced by Illumina NGS sequencing. Sequencing data were analyzed using Cell Ranger analysis, and results visualized by Loupe Cell Browser. CD4+ T-cells were clustered and classified based on Gene expression and dCODE Dextramer. Clusters of EBV- and TT- specific T-cells were detected by specific dCODE Dextramer.

B) EBV- and TT- specific T-cells were identified in both unsorted and sorted CD4+ T-cell populations.

C) The top 3 identified EBV- and TT-specific TCR clonotypes were conserved across unsorted and sorted CD4+ T-cell populations. For each clonotype, the corresponding number of cells is reported.

Conclusion

In this study we have:

- Set up a protocol to stimulate antigen-specific T cell response in a CD4+ T cell population.
- Identified antigen-specific T-cell populations in the CD4+ T cell sample.
- Identified TCR clonotypes in the antigen-specific populations, and sequenced the relative VDJ combinations, using MHC II dCODE Dextramer in combination with the Feature Barcoding technology.

For the first time, DR1-restricted TCR clonotypes have been identified in EBV- and TT-specific T cell populations.

The identification of novel MHC II-restricted TCR sequences offers the possibility to progress immunological research and immunotherapeutic development.