

Novel insights into diabetes and adaptive immunity provided by the $\alpha 4$ -integrin deficient NOD mouse

Oulghazi S. *et al.* Adaptive Immunity and Pathogenesis of Diabetes: Insights Provided by the $\alpha 4$ -Integrin Deficient NOD Mouse. *Cells* 2020, 9 (12):2597

BACKGROUND

Autoimmune type 1 diabetes (T1D) is a complex autoimmune disease characterized by the destruction of the insulin-producing β -cells in the pancreas.

In this prospective study, researchers studied the role of $\alpha 4$ (CD49d) hematopoietic integrin in T1D, using newly generated non-obese diabetic (NOD) $\alpha 4$ knock-out mice (NOD. $\alpha 4$ ^{-/-}) as a model. NOD mice exhibit spontaneous development of autoimmune T1D due to insulinitis, an inflammation of the islets of Langerhans in the pancreas.

STUDY DESCRIPTION

Goal: Evaluate the contributions of $\alpha 4$ -integrin to autoimmune diabetes using NOD. $\alpha 4$ knock-out mice and adoptive T-cell transfer experiments.

Mice from three different cohorts (pre-diabetic NOD, diabetic NOD, NOD. $\alpha 4$ ^{-/-}) were assessed for adaptive cellular and humoral immune responses against islet autoantigens and subjected to microbiota analyses. Diabetes was diagnosed based on recurrent hyperglycemia (blood glucose level >200 mg/dL).

MHC Dextramer[®] H-2 Kd/VYLKTNVFL (from islet-specific glucose-6-phosphatase catalytic subunit-related protein, IGRP) was used for detecting IGRP-autoreactive T cells by flow cytometry. For negative control, the same staining with leukocytes from MHC-disparate C57Bl/6 mice was performed.

RESULTS

- NOD. $\alpha 4$ ^{-/-} mice were completely protected from autoimmune diabetes
- NOD. $\alpha 4$ ^{-/-} mice developed islet-specific T-cells and antibodies, albeit quantitatively less than $\alpha 4$ ⁺ counterparts (**Fig.1**)
- Transplantation with isogenic $\alpha 4$ ^{-/-} bone marrow prevented progression to T1D of pre-diabetic NOD. $\alpha 4$ ⁺ mice despite significant pre-existing islet cell injury
- Transfer of $\alpha 4$ ⁺/CD3⁺, but not $\alpha 4$ ⁺/CD4⁺ splenocytes from diabetic to NOD. $\alpha 4$ ^{-/-} mice, induced diabetes with short latency
- Microbiota of diabetes-resistant NOD. $\alpha 4$ ^{-/-} and pre-diabetic NOD. $\alpha 4$ ⁺ mice were identical and are distinct from diabetic NOD mice.




Fig. 1. Diabetes and adaptive immune responses against islet antigens. MHC Dextramer[®] H-2Kd/VYLKTNVFL binding to CD8⁺ T cells was assessed by flow cytometry.

Representative dot plots for NOD, NOD. $\alpha 4$ ^{-/-} and MHC-disparate C57Bl/6 (negative control) mouse blood (left) and quantitative analysis (right), where Dextramer[®]-positive events in negative control blood were subtracted as background ($n = 9$ per group; $p < 0.05$).

CONCLUSIONS

- "NOD. $\alpha 4$ ^{-/-} mice are diabetes resistant despite developing adaptive immunity, albeit attenuated, against islet autoantigens"
- $\alpha 4$ is a promising target for primary or secondary prevention of human type 1 diabetes
- MHC I Dextramer[®] is a sensitive tool that can be used for the detection and enumeration of islet-specific autoreactive T cells in mouse model
- MHC I Dextramer[®] reagents can support and advance scientific discoveries in autoimmunity.