

CD8⁺ T-Cell Signature in Acute SARS-CoV-2 Infection Identifies Memory Precursors with dCODE Dextramer®

S. Adamo et al. *CD8⁺ T Cell Signature in Acute SARS-CoV-2 Infection Identifies Memory Precursors* [2021] *BioRxiv* [doi: <https://doi.org/10.1101/2021.07.22.453029>]

BACKGROUND

In this study, Adamo *et al.* investigate the signature of SARS-CoV-2 specific long-lived memory CD8⁺ T cells using spectral flow cytometry combined with cellular indexing of transcriptomes and T-cell receptor (TCR) sequences in PCR confirmed COVID-19 patients.

STUDY DESCRIPTION


SARS-CoV-2-specific CD8⁺ T cells were profiled directly *ex vivo* by single-cell RNA (scRNA) and TCR sequencing in PBMCs from 33 patients during acute COVID-19 as well as at six and 12 months after primary infection. SARS-CoV-2-specific CD8⁺ T cells were detected by using HLA-A*01:01, HLA-A*11:01, and HLA-A*24:02 MHC-I dCODE Dextramer® reagents.

RESULTS


scRNAseq analysis of the temporal phenotypic changes separated the SARS-CoV-2 specific CD8⁺ T cells into 12 distinct clusters based on gene expression (**Fig. 1A**). Clusters 1, 2, and 12 dominated the acute phase, while cluster 11 dominated the recovery phase (**Fig. 1B**). Clusters 1, 2, and 12 corresponded to cytotoxic, activated, and proliferating cells, while cluster 11 was enriched in IFN, TNF, and LT- α genes.

TCR sequencing revealed that the number of persistent clones decreased over time (**Fig. 1C**) and that those were enriched in genes involved in IFN- γ response, IFN- α response, and TNF signaling. Non-persisting cells were enriched in mTOR signaling genes and genes related to mitosis (**Fig. 1D**).


A) Clustering of SARS-CoV-2 specific CD8⁺ T cells



B) SARS-CoV-2 Clonotypes based on TCR sequencing



C) COV2-Dex+ CD8⁺ T cell clones (acute phase)



D) Persistent vs. non-persistent




Fig. 1: Temporal phenotypic transcriptional changes in SARS-CoV-2 CD8⁺ T cells. A) scRNASeq analysis of SARS-CoV-2 specific CD8⁺ T cells revealed 12 clusters with certain clusters dominating either acute or recovery phase. B) Clonotypes of SARS-CoV-2 specific CD8⁺ T cells determined by TCR sequencing. C) Distribution of persistent and non-persistent clones in acute phase. D) Gene expression in persistent and non-persistent SARS-CoV-2 CD8⁺ T cells.

CONCLUSIONS

- SARS-CoV-2 specific CD8⁺ T cells detected using dCODE Dextramer® reagents exhibited a T_{effector} phenotype during acute disease expressing genes related to cytotoxic, activated, and proliferating cells. During recovery, a progressive switch to a memory CD8⁺ T-cells phenotype with expression of genes related to IFN, TNF, and LT- α genes was observed in the SARS-CoV-2 specific CD8⁺ T cells
- Persisting SARS-CoV-2 CD8⁺ T-cell clones decreased over time but were enriched in genes involved in the IFN- γ response, IFN- α response, and TNF signaling