

Cancer Neopeptides for Immunotherapy: Developing Strategies for Adoptive T-Cell Therapy

[Wickström S.L. et al. Cancer neopeptides for immunotherapy: discordance between tumor-infiltrating T cell reactivity and tumor MHC peptidome display. Front. Immunol. 2019, 10:2766.](#)

Tumor-infiltrating lymphocytes (TILs) can recognize shared tumor-associated antigens (TAAs) and mutation-derived neopeptides specific to the tumor or individual patient. Therefore, adoptive T-cell therapy based on one of these features of TILs is an attractive prospect for cancer treatment.

This study highlights the challenges of activating autologous TILs by patients' neopeptides and suggests strategies for better neoantigen-targeted immunotherapy.

STUDY DESCRIPTION

Goal: to study immunogenic neopeptides using two distinct strategies: peptide reactivity of T cells isolated from the tumor and mass spectrometry detection of neopeptides presented on the surface of the tumor cells.

- Exome sequencing and epitope prediction were performed from tumor cell lines from two HLA-A2*0201 melanoma patients (KADA and ANRU) with strong tumor reactivity of TILs
- Activation of TILs was analyzed by IFN- γ ELISA or flow cytometry
- Recognition of melanoma TAAs and neopeptides by TILs was assessed using either MHC Dextramer[®] Melanoma Panel or neopeptide-specific MHC I Dextramer[®] staining by flow cytometry (**Fig.1.**).

RESULTS

- TILs recognized 5/181 and 3/49 of the predicted neopeptides
- TILs were unable to recognize the MS-defined neopeptides AGPS and ENC1 detected on tumor cells of KADA
- In KADA, TILs were detected to be specific for KDELR2, MYLIP, and SVIL epitopes, but not for WDR75
- In ANRU, all 3 predicted neopeptide-specific TILs (ETV6 9mer, ETV6, 10mer, and NUP210) resulted in well-defined populations, comparable to the staining of MART-1.




Fig.1. Detection of activated TILs was assessed by staining with the neopeptide-specific MHC I Dextramer[®] reagents. The MHC Dextramer[®] Melanoma Panel was used for staining MART-1-specific TILs by flow cytometry. FLU-specific Dextramer[®] was served as a positive control.

CONCLUSIONS

- The recognition of neopeptides is highly specific, and tolerance to wild-type antigens has not been broken since they could not activate TILs
- The ability of TILs to recognize autologous tumor cells can be efficiently assessed using either Melanoma Dextramer[®] Panel that recognize shared TAAs or neopeptide-specific MHC I Dextramer[®] reagents for mutation-derived neopeptides.