

Protocol for Preparation of MHC I-Peptide Monomer and Fluorescent U-Load dCODE Dextramer[®] MHC I reagents

Background

easYmer[®] are highly active formulations of MHC I molecules, which can be used to generate specific MHC-peptide monomers by loading your own peptide of choice. The MHC-peptide monomers can easily be dextramerized with U-Load dCODE Dextramer[®] and used for profiling and quantification of antigen-specific T cells in cell sample. The U-Load dCODE Dextramer[®] reagent is a DNA barcode labelled Dextramer backbone with a unique barcode label for each specificity. U-Load dCODE Dextramer[®] reagent comes with labels applicable for different applications:

- U-Load dCODE Dextramer[®] (RiO) for multi-omics analysis of antigen-specific T cells, compatible with BD's Rhapsody single cell analysis platform
- U-Load dCODE Dextramer[®] (10x) – for multi-omics analysis of antigen-specific T cells, compatible with 10x Chromium Single Cell Gene Expression system
- U-Load dCODE Dextramer[®] (HiT) - for epitope discovery, neoantigen screening, designed for multiplexing using PCR and NGS,

The easYmer[®] and U-Load Dextramer[®] technologies are highly flexible and suitable for screening of single epitopes in many samples or screening of large numbers of different epitopes in parallel.

The easYmer[®] technology also allow evaluation of peptide binding to MHC I by assaying proper refolding of peptide loaded complexes.

Materials required

easYmer[®]
 U-Load dCODE Dextramer[®] (RiO, 10x, HiT)
 Peptide
 DMSO (e.g., Sigma cat# D2650)
 ddH₂O

Preparation of MHC-peptide monomer

1. Dilute Peptide (easYmer[®] Control Peptide or peptide of choice) to 100 µM in ddH₂O. Keep on ice from this step.
2. To prepare easYmer[®] MHC-peptide monomer for dextramerization, mix the reagents in Table A according to the listed sequence in a 1.5 mL tube or 96-well U-bottom plate. The listed amounts will be enough to make 10, 20, or 50 tests of U-Load dCODE Dextramer[®] MHC I reagent.
3. *Optional: To evaluate folding of peptide with MHC make a smaller volume of a positive and negative control, i.e., easYmer[®] loaded with the included control peptide and no peptide as listed in Table A.*

Table A

Reagents	10 tests	20 tests	50 tests	Positive Control	Negative control
ddH ₂ O	3 µL	6 µL	15 µL	2.5 µL	3 µL
Peptide (100 µM)	2 µL	4 µL	10 µL	0.5 µL	-
easYmer [®] Loading Buffer	5 µL	10 µL	25 µL	3 µL	3 µL
easYmer [®] (3 µM)	20 µL	40 µL	100 µL	3 µL	3 µL

For research use only. Not for use in diagnostic or therapeutic procedures.

Total volume	30 µL	60 µL	150 µL	9 µL	9 µL
--------------	-------	-------	--------	------	------

4. Mix by pipetting gently – *be careful not to form bubbles.*
5. Make a short clearing spin to collect all materials in the bottom of the tube and incubate at 18°C for 48 h.
6. Make a short clearing spin to collect all material in the bottom of the tube. 2000 nM folded MHC-peptide monomer is now ready for further processing or placed at -20°C for long-term storage.
7. Proceed to step 8 to evaluate peptide binding or continue to dextramerization starting at step 9.
To determine the efficiency of the easYmer® MHC-peptide folding take 2 µL of the prepared MHC-peptide monomer (2000 nM) and dilute to 500 nM by adding 6 µL of Loading buffer. Then follow the “Protocol for Flow cytometry-based assay of peptide- easYmer® MHC I complex formation” (page 3).

Preparation of U-Load dCODE Dextramer® MHC I

8. To dextramerize the peptide-loaded easYmer® MHC monomers (MHC-peptide monomer) with U-Load dCODE Dextramer® reagent, mix the reagents in Table B in a 1.5 mL tube:

Table B

Reagents	10 tests	20 tests	50 tests
MHC-peptide monomer (2 µM)	27 µL	54 µL	135 µL
U-Load dCODE Dextramer®	12 µL	24 µL	60 µL
<i>incubate for 30 min at RT in the dark</i>			
U-Load dCODE Dextramer® Dilution Buffer	11 µL	22 µL	55 µL
Total volume U-Load dCODE Dextramer® MHC I	50 µL	100 µL	250 µL

9. Store the fluorescent DNA barcode labelled U-Load dCODE Dextramer® MHC I reagents at 2-8°C in the dark until use. *5 µL per test is recommended for staining of 1-10 x 10⁶ PBMC.*

Analysis of antigen-specific T cells

Antigen-specific T cells can be analyzed on three sequencing platforms:

- U-Load dCODE Dextramer® (HiT) for epitope discovery, neoantigen screening and immune monitoring follow the protocol “dCODE Dextramer® (HiT) Package Insert”
- U-Load dCODE Dextramer® (10x) for multi-omics analysis of antigen-specific T cells and TCR sequencing using 10x Chromium Single Cell Gene Expression follow the protocol “dCODE Dextramer® (10x compatible) Staining Protocol”
- U-Load dCODE Dextramer® (RiO) for multi-omics analysis of antigen-specific T cells and gene expression profiling using the BD Rhapsody™ Single-Cell Analysis System follow the protocol “dCODE Dextramer® (RiO) Package Insert” (<https://www.immudex.com/resources/protocols/>)


For research use only. Not for use in diagnostic or therapeutic procedures.

Flow Cytometry-Based Assay of Peptide-MHC I Complex Formation

Background After completed peptide loading of easYmer[®] monomer (step 8 in the protocol “Preparation of MHC I-peptide monomer and fluorescent U-Load Dextramer[®] MHC I”), proper folding of MHC-peptide monomers can be evaluated using the assay described in this procedure.

Procedure: evaluation of peptide-MHC I complex formation

1. The MHC-peptide monomers to be assayed are prepared according to protocol: “Preparation of MHC I-peptide monomer and fluorescent U-Load Dextramer[®] MHC I”. The peptide loading setup should include a positive control (provided with the easYmer[®]): a peptide that is known to support folding of the MHC molecule of interest, and a negative control: without peptide.
2. Prepare sufficient dilution buffer (PBS, 5% glycerol) for the whole assay.
3. Dilute each of the MHC-peptide monomers to give 75 μ L of a 40 nM solution (e.g., for a 500 nM complex: 6 μ L folded complex in 69 μ L dilution buffer).
4. For all samples and positive and negative controls, transfer 50 μ L of this pre-dilution (prepared in step 3.) to the first tube. Make three subsequent serial 3-fold dilutions (50 μ L in 100 μ L dilution buffer), according to the figure below.



5. Transfer 40 μ L of each of these dilutions to the wells in a U-bottom shape 96-well plate, as suggested below. Also prepare a background well (BLANK): transfer 40 μ L of dilution buffer
6. Prepare a sufficient volume of a 45-fold dilution of the streptavidin-coated beads (6-8 μ m; SVP-60-5) in dilution buffer. Transfer 20 μ L of the diluted bead suspension to each well.

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLANK		P-1		S1-1		S3-1		S5-1		S7-1	
B			P-2		S1-2		S3-2		S5-2		S7-2	
C			P-3		S1-3		S3-3		S5-3		S7-3	
D												
E			N-1		S2-1		S4-1		S6-1		S8-1	
F			N-2		S2-2		S4-2		S6-2		S8-2	
G			N-3		S2-3		S4-3		S6-3		S8-3	
H												

BLANK : No complex


P1-3 : Positive control dilutions (HLA with know peptide)

N1-3 : Negative control dilutions (HLA without peptide)

S1-S8 : Sample dilutions (complexes to evaluate)

7. Mix well and seal the plates with Sealing Tape to avoid well to well contamination.
8. Incubate the plate on a rocking table at 37°C for 1h.
9. Remove the Sealing Tape and wash by adding 160 µL FACS buffer.
10. Spin the plate at 700xg for 3 min and flip out the supernatant.
11. Resuspend the beads in 200 µL FACS buffer.
12. Spin the plate at 700xg for 3 min and flip out the supernatant.
13. Wash two more times by repeating step 11. and 12.
14. During the above washing steps, prepare a 200-fold dilution of the PE-labeled anti-human β2m monoclonal antibody BBM.1 in FACS buffer.
15. Resuspend the beads in 50 µL antibody solution per well.
16. Incubate the plate for 30 min at 4°C.
Wash by adding 150 µL FACS buffer. Spin the plate at 700g for 3 min and flip out the supernatant.
17. Resuspend the beads in 200 µL FACS buffer. Spin the plate at 700xg for 3 min and flip out the supernatant.
18. Wash two more times by repeating step 17. and 18.
19. Resuspend the beads in 200 µL FACS buffer and analyze on a Flow cytometer.

Example of the Flow cytometry-based assay:



Flow cytometry-based detection of 4 different peptide-HLA-A*02:01 monomers. MHC-peptide monomers of A*02:01 and 4 different peptides, and a negative control (No Peptide), were folded. CMV pp65 495-503 (NLVPMVATV) a known HLA-A*02:01 restricted epitope was used as positive control. The three other peptides are based on their A*02:01 binding stability categorized as good binder ($T_{1/2}$ 6.5h), intermediate binder ($T_{1/2}$ 3.5h), and low binder ($T_{1/2}$ 0.7h). Three dilutions of the folded complexes were analysed in the flow cytometry-based assay. The X-axis gives the complex concentration if complete folding is achieved.

For research use only. Not for use in diagnostic or therapeutic procedures.