

## MHC Dextramer® Staining Protocol – Single Reagent on PBMCs

**Intended use** Staining of antigen-specific T cells using a single fluorochrome-labelled MHC Dextramer® reagent.

**Materials Provided** MHC Dextramer® PE, APC or FITC

**Materials Required (not provided)** Stain and wash buffer: PBS, 1-5% FCS, pH 7.4  
Antibodies identifying relevant surface markers (e.g., CD3, CD4, CD8 and optionally other desired antibodies). The optimal choice of antibodies and fluorochromes depends on the flow cytometer and experimental setup.

**Procedure**

1. Thaw and prepare PBMC sample and resuspend  $1-3 \times 10^6$  PBMC in 50  $\mu$ L stain and wash buffer.
2. Centrifuge MHC Dextramer® reagents at 10,000 x g for 1 min.
3. Add 10  $\mu$ L MHC Dextramer® to the cell sample and vortex briefly.
4. Incubate at room temperature in the dark, for 10 min. *MHC II Dextramer® reagents require longer incubation. See Procedural notes.*
5. Add relevant antibodies in the volume/concentration recommended by the provider:
  - a. If staining with MHC I Dextramer® reagents, use anti-CD8, anti-CD3, and optionally other phenotype markers.
  - b. If staining with MHC II Dextramer® reagents, use anti-CD4, anti-CD3, and optionally other phenotype markers.
6. Incubate at room temperature in the dark for 20 min.
7. Wash cells by adding 2 mL stain and wash buffer. Centrifuge at 300 x g for 5 min. and remove the supernatant. Repeat washing for a total of 2 washes.
8. Resuspend the pellet in 300  $\mu$ L or desired volume stain and wash buffer suitable for your flow cytometer.
9. Proceed to analyze the samples on a flow cytometer or store at 2-8°C in the dark. For optimal results, do not store the samples longer than 2 hours before acquisition. *Alternatively, samples can be fixated. For fixated samples, see Procedural notes.*

**Procedural notes**

1. If you intend to stain with multiple MHC Dextramer® reagents, follow the protocol "MHC Dextramer® Staining Protocol – Multiple Reagents". (<https://www.immudex.com/resources/protocols/>)
2. MHC Dextramer® staining can be performed on whole blood.
3. Protocol step 4: When staining with MHC II Dextramer® reagents, a mixed MHC I and II Dextramer® pool or identifying low affinity targets, stain for 30 min.
4. Protocol step 7: Staining can be performed using 96-well microtiter plates. In that case after antibody incubation make 4 sequential washes using 200  $\mu$ L stain and wash buffer per well. Centrifuge at 300 x g for 5 min between each wash and remove supernatant.

5. Protocol step 9: MHC Dextramer<sup>®</sup> stained cells can be fixated using 2% formaldehyde in PBS. Fixated samples can be stored in the dark for up to 24 hours before acquisition on a flow cytometer. Fixated samples should be washed and resuspended in stain and wash buffer prior to acquisition on a flow cytometer.
6. Always keep MHC Dextramer<sup>®</sup> reagents stored at 2-8°C in the dark – the plastic vial only partially protects the reagents against light.

**Storage** Store in the dark at 2-8°C.

**Precautions** Contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. For professional users.

**Symbols** See [www.immudex.com/symbols](http://www.immudex.com/symbols) for explanation of symbols.

**Technical support** E-mail: [customer@immudex.com](mailto:customer@immudex.com)  
Telephone: +45 3110 9292 (Denmark)

**Manufacturer** Immudex, Bredevej 2A, DK-2830 Virum, Denmark