

## Protocol for preparation and loading of easYmers® MHC I-peptide monomer onto U-Load dCODE Dextramer®

### Background

easYmers® powered by immunAware is a formulation of peptide-receptive MHC I monomer, which can be used to generate specific MHC I-peptide monomer by loading your peptide of choice. The easYmers® MHC I-peptide monomer can easily be loaded onto U-Load dCODE Dextramer® and used for characterization and quantification of antigen-specific T cells in a cell sample by next-generation sequencing (NGS) or single-cell multi-omics. U-Load dCODE Dextramer® is a DNA barcode labelled Dextramer® with a unique DNA barcode for each specificity. In addition, U-Load dCODE Dextramer® is labeled with PE for cell-sorting purposes. U-Load dCODE Dextramer® comes with DNA barcodes applicable for different applications:

- U-Load dCODE Dextramer® (HiT) - for epitope discovery, neoantigen screening, designed for multiplexing using PCR and NGS
- U-Load dCODE Dextramer® (RiO) for the detection of antigen-specific CD8<sup>+</sup> or CD4<sup>+</sup> T cells with additional information of gene expression, surface marker expression, and TCR sequence by single-cell multi-omics using the BD Rhapsody™ Single-Cell Analysis System
- U-Load dCODE Dextramer® (10x) for the detection of antigen-specific CD8<sup>+</sup> and CD4<sup>+</sup> T cells with additional information on gene expression, surface marker expression, and V(D)J sequences using the 10x Chromium™ platforms

The easYmers® and U-Load dCODE Dextramer® technologies are highly flexible and suitable for screening of single epitopes in many samples or screening of large numbers of different epitopes in parallel. The easYmers® technology also allows evaluation of peptide binding to MHC I by assaying proper refolding of peptide-loaded monomer.

### Materials required

The materials listed here are required for preparation of easYmers® MHC-peptide monomer and U-Load Dextramer® MHC I.

easYmers®  
easYmers® loading buffer  
easYmers® positive control peptide  
U-Load dCODE Dextramer® (HiT, RiO, 10x)  
U-Load dCODE Dextramer® dilution buffer

### Materials required (not provided)

The materials listed here are required for preparation of easYmers® MHC I-peptide monomer and U-Load dCODE Dextramer® MHC I and for the flow cytometry-based assay for evaluation of proper folding of easYmers® MHC I-peptide monomer.

Peptide of choice  
DMSO (e.g., Sigma cat.# D2650)

Dilution buffer (PBS, 5% glycerol)  
 FACS buffer (PBS, 1% BSA (or FCS), 0.01% NaN<sub>3</sub>)  
 Streptavidin-coated beads (Spherotech cat.# SVP-60-5)  
 Anti-human  $\beta_2m$  BBM.1-PE (Santa Cruz cat.# sc-13565 PE)

### Protocol steps and timing

Experimental workflow using the easYmers<sup>®</sup> and U-Load dCODE Dextramer<sup>®</sup> and estimated time to complete each step.



### I. Preparation of easYmers<sup>®</sup> MHC I-peptide monomer

1. Reconstitute your peptides of interest according to the manufacturer's instructions.
2. Dilute Peptide (easYmers<sup>®</sup> control peptide or peptide of interest) to 100  $\mu$ M in ddH<sub>2</sub>O. Keep on ice from this step on.
3. To prepare easYmers<sup>®</sup> MHC I-peptide monomer, mix the reagents in Table A according to the listed sequence in a 1.5 mL tube or 96-well U-bottom plate. The listed amounts will be enough to make 10, 20, or 50 tests of U-Load dCODE Dextramer<sup>®</sup> MHC I.

*Optional: To evaluate the peptide loading efficiency make a smaller volume of the easYmers<sup>®</sup> positive and the negative control (no peptide), i.e., easYmers<sup>®</sup> loaded with the included easYmers<sup>®</sup> positive control peptide or no peptide as listed in Table A.*

**Table A**

Reagents	10 tests	20 tests	50 tests	Positive Control	Negative control
ddH <sub>2</sub> O	3 $\mu$ L	6 $\mu$ L	15 $\mu$ L	2.5 $\mu$ L	3 $\mu$ L
Peptide (100 $\mu$ M)	2 $\mu$ L	4 $\mu$ L	10 $\mu$ L	0.5 $\mu$ L	-
easYmers <sup>®</sup> Loading Buffer	5 $\mu$ L	10 $\mu$ L	25 $\mu$ L	3 $\mu$ L	3 $\mu$ L
easYmers <sup>®</sup> (3 $\mu$ M)	20 $\mu$ L	40 $\mu$ L	100 $\mu$ L	3 $\mu$ L	3 $\mu$ L
<b>Total volume of easYmers<sup>®</sup> MHC I-peptide monomer (2 <math>\mu</math>M)</b>	<b>30 <math>\mu</math>L</b>	<b>60 <math>\mu</math>L</b>	<b>150 <math>\mu</math>L</b>	<b>9 <math>\mu</math>L</b>	<b>9 <math>\mu</math>L</b>

4. Mix by pipetting gently – *be careful not to form bubbles.*
5. Briefly centrifuge to collect all materials in the bottom of the tube and incubate at 18°C for 48 h.
6. Briefly centrifuge to collect all material in the bottom of the tube. 2  $\mu$ M folded MHC I-peptide monomer are now ready for loading onto U-Load dCODE Dextramer<sup>®</sup> backbone or can be stored at -20°C for long-term storage.

7. Proceed to page 4 to evaluate peptide loading efficiency or continue to load onto U-Load dCODE Dextramer®.

## II. Loading of U-Load dCODE Dextramer® MHC I

8. To load the easYmers® MHC I-peptide monomer onto U-Load dCODE Dextramer®, mix the reagents in Table B in a 1.5 mL tube:

**Table B**

Reagents	10 tests	20 tests	50 tests
easYmers® MHC I-peptide monomer (2 µM)	27 µL	54 µL	135 µL
U-Load dCODE Dextramer® (PE)	12 µL	24 µL	60 µL
<i>incubate for 30 min at RT in the dark</i>			
U-Load dCODE Dextramer® Dilution Buffer	11 µL	22 µL	55 µL
<b>Total volume</b> <b>U-Load dCODE Dextramer® MHC I</b>	<b>50 µL</b>	<b>100 µL</b>	<b>250 µL</b>

9. Store the fluorescent U-Load dCODE Dextramer® MHC I reagents at 2-8°C in the dark until use.

## III. Staining Procedures & Sequencing Workflows

- For U-Load dCODE Dextramer® (HiT): See [www.immudex.com/Protocols/HiT](http://www.immudex.com/Protocols/HiT)  
 For U-Load dCODE Dextramer® (RiO): See [www.immudex.com/Protocols/RiO](http://www.immudex.com/Protocols/RiO)  
 For U-Load dCODE Dextramer® (10x): See [www.immudex.com/Protocols/10x](http://www.immudex.com/Protocols/10x)


## Optional: Flow cytometry-based quality control assay for determination of peptide loading efficiency

### Background

After easYmers<sup>®</sup> MHC I-peptide monomerization (step 6 in the protocol), the relative peptide-loading efficiency can be determined by comparing your peptide of interest to the negative and positive loading controls using this assay. The negative loading control is empty easYmers<sup>®</sup> (no peptide). The positive loading control peptide is specific to and provided with the easYmers<sup>®</sup> you purchase. If this is your first time testing a particular easYmers<sup>®</sup> MHC I-peptide combination, this assay is highly recommended.

### Procedure: evaluation of easYmers<sup>®</sup> MHC I-peptide monomer formation

1. Prepare a sufficient volume of dilution buffer (PBS, 5% glycerol).
2. To determine the efficiency of the easYmers<sup>®</sup> MHC I-peptide folding take 3  $\mu\text{L}$  of the prepared easYmers<sup>®</sup> MHC I-peptide monomer (1  $\mu\text{M}$ ) and dilute to 500 nM by adding 3  $\mu\text{L}$  of dilution buffer.
3. Dilute each of the easYmers<sup>®</sup> MHC I-peptide monomer to give 75  $\mu\text{L}$  of a 40 nM solution (e.g., for a 500 nM monomer: 6  $\mu\text{L}$  folded monomer in 69  $\mu\text{L}$  dilution buffer).
4. For all samples and positive and negative loading controls, transfer 50  $\mu\text{L}$  of this pre-dilution (prepared in step 3) to the first tube. Make three subsequent serial 3-fold dilutions (50  $\mu\text{L}$  in 100  $\mu\text{L}$  dilution buffer), according to the figure below.



5. Transfer 40  $\mu\text{L}$  of each of these dilutions to the wells in a U-bottom shape 96-well plate, as suggested below. Also, prepare a background well (BLANK): 40  $\mu\text{L}$  of dilution buffer (no beads or antibody will be added to this well).

6. Prepare a sufficient volume of a 45-fold dilution of the streptavidin-coated beads in dilution buffer. Transfer 20  $\mu$ L of the diluted bead suspension to each well.

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLANK		P-1		S1-1		S3-1		S5-1		S7-1	
B			P-2		S1-2		S3-2		S5-2		S7-2	
C			P-3		S1-3		S3-3		S5-3		S7-3	
D												
E			N-1		S2-1		S4-1		S6-1		S8-1	
F			N-2		S2-2		S4-2		S6-2		S8-2	
G			N-3		S2-3		S4-3		S6-3		S8-3	
H												

**BLANK** : No complex


**P1-3** : Positive control dilutions (HLA with know peptide)

**N1-3** : Negative control dilutions (HLA without peptide)

**S1-S8** : Sample dilutions (complexes to evaluate)

7. Mix well and seal the plates with sealing tape to avoid well to well contamination.
8. Incubate the plate on a rocking table at 37°C for 1 h.
9. Remove the sealing tape and wash by adding 160  $\mu$ L FACS buffer.
10. Spin the plate at 700 x g for 3 min and discard the supernatant.
11. Resuspend the beads in 200  $\mu$ L FACS buffer.
12. Spin the plate at 700 x g for 3 min and discard the supernatant.
13. Wash two more times by repeating step 10 and 12.
14. During the above washing steps, prepare a 200-fold dilution of the PE-labeled anti-human  $\beta$ 2m monoclonal antibody BBM.1 in FACS buffer.
15. Resuspend the beads in 50  $\mu$ L antibody solution per well.
16. Incubate the plate for 30 min at 4°C.
17. Wash by adding 150  $\mu$ L FACS buffer. Spin the plate at 700 x g for 3 min and discard the supernatant.
18. Resuspend the beads in 200  $\mu$ L FACS buffer. Spin the plate at 700 x g for 3 min and discard the supernatant.
19. Wash two more times by repeating step 17 and 18.
20. Resuspend the beads in 200  $\mu$ L FACS buffer and analyze on a Flow cytometer.

**Example of the Flow cytometry-based assay:**



**Flow cytometry-based detection of 4 different peptide-HLA-A\*02:01 monomers.** MHC I-peptide monomer of A\*02:01 and 4 different peptides, and a negative control (no peptide), were folded. CMV pp65 495-503 (NLVPMVATV) a known HLA-A\*02:01 restricted epitope was used as positive control. The three other peptides are based on their A\*02:01 binding stability categorized as good binder ( $T_{1/2}$  6.5h), intermediate binder ( $T_{1/2}$  3.5h), and low binder ( $T_{1/2}$  0.7h). Three dilutions of the folded monomer were analysed in the flow cytometry-based assay. The X-axis gives the monomer concentration if complete folding is achieved.