

Advancements in single cell multiomic profiling of antigen-specific T cells with dCODE Dextramer® (RiO) and BD® AbSeq Reagents on the BD Rhapsody™ Single-Cell Analysis System



Cynthia Sakofsky,¹ Kivin Jacobsen,² Vadir Lopez-Salmeron,⁴ Liselotte Brix,² Aruna Ayer¹.


¹BD Biosciences, 2350 Qume Drive, San Jose, CA 95131, ²Immudex, Fruebjergvej 3 DK-2100 Copenhagen, Denmark, ³BD Biosciences, Tullastrasse 8-12 69126 Heidelberg, Germany

Abstract

The advancement of immunotherapies and drug development relies on the characterization of antigen-specific T cells. However, due to their low frequency among T cells and their inherent weak affinity to bind known MHC-peptide complexes, detection of these cells has been difficult. Furthermore, pairing this information with the corresponding variable (V), diversity (D), and joining (J) [V(D)J] sequences of the antigen-specific T cell receptors has also been challenging due to hurdles in sequencing these large regions, particularly at the single-cell level. Here, we have expanded our previous work of combining two powerful technologies, Immudex® dCODE Dextramer® (RiO) Reagents and the BD Rhapsody™ Single-Cell Analysis System, to detect and characterize low-frequency antigen-specific T cells, including the full sequences of the V(D)J gene segments of the T cell receptors, as well as profile transcriptome and protein expression. Specifically, thousands of sorted PBMCs were multiplexed to provide high-throughput detection of individual antigen-specific CD8⁺ T cells in combination with the corresponding full V(D)J sequences of the T cell receptors. In addition, we simultaneously obtained gene expression data for over 350 immune-related mRNAs as well as cell phenotypes using a panel of 15 cell surface BD® AbSeq Protein Markers. Together these data can be used to define T cell phenotypes associated with T cell activation states, alongside antigen specificity of enriched CD8⁺ dCODE Dextramer® positive antigen-specific T cells from PBMC populations. This study showcases the importance of multiomic analysis for high-resolution T cell profiling that has broader implications and utility in immuno-oncology, infectious diseases and autoimmunity.

Methods

Seamless integration of dCODE Dextramer® (RiO) technology into the BD® Rhapsody AbSeq + SMK + TCR Full Length Workflow



dCODE Dextramer® reagents were used alongside BD® AbSeqs to sensitively identify antigen specific T cells and profile cell surface proteins




Figure 2. dCODE (RiO) reagents and BD® AbSeqs were used simultaneously to identify antigen-specific T cells in hPBMCs and profile cell surface proteins. (A) tSNE plot of an experiment using ~ 5K sorted dCODE Dextramer positive CD8⁺ cells with cell-type distribution shown. (A) tSNE plots highlighting the expression of a 13-plex dCODE Dextramer panel including 8 dCODE reagents with epitopes and 5 negative controls. (B) tSNE plots showing the expression of a 15-plex BD AbSeq panel that range in expression level. (C) Full length translated TCR sequence of the top six clonotypes.

Full length TCR sequences were identified in antigen specific CD8⁺ T cells and were linked to specific antigen epitopes using dCODE® (RiO) reagents




Figure 3. Specific antigen epitopes linked to full length TCR sequences. (A) The six most frequent TCR CDR3 clonotypes were identified among CD8⁺ dCODE Dextramer positive antigen-specific T cells along with their mRNA gene expression profiles of 26 high expressing genes from the BD Rhapsody™ Immune Response panel (397 genes). (B) Differential gene expression plots of the 13 dCODE Dextramer panel for each of the top six clonotypes show a high association to a single dCODE Dextramer. (C) Full length translated TCR sequence of the top six clonotypes.

BD Rhapsody sample multiplexing is compatible with dCODE Dextramer® (RiO) reagents and provides a solution for high throughput detection of antigen-specific T cells




Figure 4. Sample multiplexing with BD Rhapsody Sample Multiplexing Kit is compatible with dCODE Dextramer® technology. Three cell samples stained with dCODE (RiO) reagents and BD AbSeq were multiplexed with BD SMK and loaded into a single BD Rhapsody cartridge. Two of the samples were CD8⁺ dCODE⁺ sorted cells derived from the same stock but independently stained with a 15-plex AbSeq and 13-plex dCODE (RiO) panel to evaluate reproducibility of AbSeq and dCODE performance as shown in (4B) Note only antigen dCODEs are shown, no negative controls. The third cell sample was an enriched CD4⁺ cell population that was stimulated with EBV and tetanus toxoid (TT) peptides and stained with a 15-plex AbSeq and 3-plex Dextramer panel. AbSeq CD4 expression and antigen-specific TCR expression detected by dCODE (RiO) reagents are shown in (4C).

Conclusions

- * dCODE Dextramer® (RiO) reagents can be seamlessly integrated into the BD AbSeq, Sample Tag, Targeted mRNA and TCR Full Length workflow to profile hPBMCs.
- * Antigen-specific T cells along with their corresponding full length TCR sequences were identified by viral specific HLA-peptide complexes displayed on dCODE Dextramer® (RiO) reagents used in combination with BD TCR Full Length assays.
- * Sample multiplexing of 3 different samples showed reproducibility of AbSeq and dCODE (RiO) performance as well as sensitive detection of dCODE (RiO) reagents from different cell types (CD4⁺ and CD8⁺ cells).

Class 1 Laser Product. For Research Use Only. Not for use in diagnostic or therapeutic procedures. Trademarks are the property of their respective owners. © 2016 BD, BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company. 23-XXXX-XX