

MHC Dextramer® Staining Protocol

Products

Cat No. Wxxxxxx/Fxxxxxx/XDxxxxx
 MHC I Dextramer® FITC, PE, or APC, cat# WBxxxxx / JDxxxxx
 MHC II Dextramer® FITC, PE, or APC, cat# FBxxxxx
 CD1d Dextramer® FITC, PE, or APC, cat# XDxxxxx / YDxxxxx
 MR1 Dextramer® FITC, PE, or APC, cat# ZAxxxxx
 Collectively denominated as Dextramer®

Recommended use

Staining of antigen-specific T cells, NKT or MAIT cells using one or more fluorochrome-labelled MHC Dextramer® reagents in one sample.

Materials Provided

MHC Dextramer® PE, APC and/or FITC
 And/or CD1d Dextramer® PE, APC, and/or FITC
 And/or MR1 Dextramer® PE, APC, and/or FITC
 Collectively denominated as Dextramer®

Materials Required (not provided)

4 mL Falcon disposable 12 x 75-mm test tubes or equivalent
 LoBind® Eppendorf tubes or equivalent
 Stain and wash buffer: PBS, 1-5% FCS, pH 7.4
 100 µM d-Biotin (e.g. Avidity, cat# BIO200) diluted in PBS, pH 7.4
 10x PBS, pH 7.4
 Antibodies identifying relevant cell surface markers:
 For CD8⁺ T, CD4⁺ T and NKT cells (e.g., CD3, CD4 and CD8).
 For MAIT cells (e.g. CD3, CD4, CD8 and CD161).
 Optionally other desired antibodies and live-dead dye^A.
 See the FAQ on immudex.com regarding [recommended antibody clones](#).
 The optimal choice of fluorochromes depends on the flow cytometer and experimental setup.

Procedure

1. Thaw and prepare PBMCs^B and resuspend 1-3 x 10⁶ PBMCs (for clonal cells, use 2-5 x 10⁴ instead) in 50 µL stain and wash buffer.
2. To prepare a pool of multiple MHC Dextramer® reagents (*calculation example can be found in Appendix 1*), mix the following reagents in an empty 1.5 mL LoBind® Eppendorf tube^C:
 - a. Add 0.2 µL of 100 µM d-Biotin^D per Dextramer® reagent.
 - b. Add 10 µL of each Dextramer® reagent.
 - c. Add 0.6 µL of 10x PBS^D per Dextramer® reagent.

NB: When staining with a single Dextramer® reagent, a and c can be omitted.
3. Vortex the Dextramer® pool briefly. The Dextramer® pool must be used directly after preparation and cannot be stored.
4. Centrifuge the pool at 10.000 x g for 1 min. to avoid transferring any potential precipitate.
5. Add the Dextramer® pool to the cell sample and vortex briefly.
6. Incubate in the dark at room temperature:
 - a. MHC I, MR1 or CD1d Dextramer® pool: 10 min. incubation^E.
 - b. MHC II Dextramer® pool: 30 min. incubation^E.

- c. Dextramer® pool comprised of a. and b.: 30 min. incubation^E.
7. Add relevant antibodies in the volume/concentration according to manufacturer's instructions:
 - a. If staining with MHC I Dextramer® reagents, use anti-CD8, anti-CD3, and optionally other phenotype markers.
 - b. If staining with MHC II Dextramer® reagents, use anti-CD4, anti-CD3, and optionally other phenotype markers.
 - c. If staining with MR1 Dextramer® reagents, use anti-CD3, anti-CD8, anti-CD4, anti-CD161 and optionally other phenotype markers.
 - d. If staining with CD1d Dextramer® reagents, use anti-CD3, anti-CD8 and anti-CD4 and optionally other phenotype markers.
8. Incubate at room temperature in the dark for 20 min.
9. Wash cells by adding 2 mL stain and wash buffer. Centrifuge at 300 x g for 5 min. and remove the supernatant. Repeat washing for a total of 2 washes^F.
10. Resuspend the pellet in desired volume of stain and wash buffer suitable for your flow cytometer.
11. Proceed to analyze the samples on a flow cytometer or store at 2-8 °C in the dark. For optimal results, do not store the samples longer than 2 hours before acquisition. Alternatively, fixed cells^G can be stored at 2-8C in dark for up to 24 hours.

Procedural notes

- A. Live-dead staining can be performed at the beginning or end of staining procedure according to manufacturer's instructions.
- B. Dextramer® staining can be performed on any cell suspensions, cell lines, TILs, or whole blood, if the cells are non-fixed. For whole-blood samples, stain with Dextramer® reagents before Red Blood Cell (RBC) lysis or use non-fixable RBC lysing solution.
- C. Always keep Dextramer® reagents stored at 2-8 °C in the dark – the plastic vial only partially protects the reagents against light.
- D. d-biotin is required to avoid artefacts in the staining. 10x PBS will balance the salt concentration of the pool.
- E. Incubation time may be increased when using a high number of reagents in pool staining and requires optimization.
- F. Staining can be performed using 96-well microtiter plates. In that case after antibody incubation make 4 sequential washes using 200 µL stain and wash buffer per well. Centrifuge at 300 x g for 5 min. between each wash and remove supernatant.
- G. Dextramer® stained cells can be fixed using 2% Methanol free formalin in PBS. Fixed samples may be washed and resuspended in stain and wash buffer prior to acquisition on a flow cytometer.

Technical support

For additional Tips & Tricks, FAQs and protocols, please visit <https://www.immudex.com/resources/> or contact our support team at customer@immudex.com
 Telephone: +45 3110 9292 (Denmark)

**Analysis
Guidelines**




Fig. 1: Flow cytometry gating strategy using MHC I Dextramer[®] to identify antigen specific T-cells from samples of thawed hPBMCs. (A-F) gating of CD8⁺ antigen specific T cells. (A) Lymphocytes were identified based on the forward (FSC) - and side scatter (SSC) profiles. **(B)** Next, doublets were excluded by gating the single cells in a side scatter height (SSC-H) & side scatter area (SSC-A) profile plot. **(C)** Dead cells were excluded according to the live-dead stain (FVS780), and the live cells were gated for further characterization. **(D)** To exclude CD4⁺ T cells and Natural killer cells (NK) (positive for CD8 but not CD3), the CD3⁺/CD4⁻ population was gated. **(E)** The CD3⁺/CD8⁺ T cells were then gated, and **(F)** subsequently, the antigen-specific population of cells were determined by comparing the results of gating the MHC I Dextramer[®] labeled or MHC I Dextramer[®] Negative Control labeled cells. **(G)** Flow cytometry plots showing CD4⁺ T helper cells labeled with MHC II Dextramer[®] or Negative Control MHC II Dextramer[®].

Appendix 1 Calculation Examples

Preparation of pools of MHC Dextramer[®] reagents for staining 1 sample:

| Examples | 100 µM d-Biotin | Total MHC Dextramer [®] Reagents | 10x PBS | Total Volume |
|--|-----------------|---|---------|--------------|
| Per each MHC Dextramer [®] | 0.2 µL | 10 µL per MHC Dextramer [®] | 0.6 µL | 10.8 µL |
| 2 MHC Dextramer [®] reagents | 0.4 µL | 20 µL MHC Dextramer [®] | 1.2 µL | 21.6 µL |
| 3 MHC Dextramer [®] reagents | 0.6 µL | 30 µL MHC Dextramer [®] | 1.8 µL | 32.4 µL |
| 10 MHC Dextramer [®] reagents | 2 µL | 100 µL MHC Dextramer [®] | 6 µL | 108 µL |

Preparation of pools of MHC Dextramer[®] reagents for staining 2 samples:

Note: When preparing a pool for more than 1 sample, we recommend preparing 20% overage of the pool, which is included in the examples below.

| Examples | 100 µM d-Biotin | Total MHC Dextramer [®] Reagents | 10x PBS | Total Volume |
|--|-----------------|---|---------|--------------|
| Per each MHC Dextramer [®] | 0.2 µL | 12 µL per MHC Dextramer [®] | 0.7 µL | 12.9 µL |
| 2 MHC Dextramer [®] reagents | 0.5 µL | 24 µL MHC Dextramer [®] | 1.4 µL | 25.9 µL |
| 3 MHC Dextramer [®] reagents | 0.7 µL | 36 µL MHC Dextramer [®] | 2.2 µL | 38.9 µL |
| 10 MHC Dextramer [®] reagents | 2.4 µL | 120 µL MHC Dextramer [®] | 7.2 µL | 129.6 µL |