

Profiling of the antigen-specific immune response using Dextramer® and dCODE® reagents in combination with flow cytometry and 10x Chromium Single-Cell Analysis System


Dilek Inekci¹, Bjarke Endel Hansen¹ and Liselotte Brix¹
¹Immudex, Copenhagen, Denmark

Introduction


Understanding the antigen-specific B and T cell responses is key for development of vaccines and targeted therapies for cancers, encompassing various stages from target discovery to monitoring the treatment efficacy to patient stratification. The Dextramer® and Klickmer® reagents allow simultaneous detection of low-frequency ag-specific B and T cells in the same workflow. For a deeper investigation, the dCODE Dextramer® and dCODE Klickmer® reagents can be used in combination with single-cell RNA sequencing, which provides a deep dive into the ag-specific B and T cells at the individual cell level giving access to BCR/TCR sequences for specific targets.

Here we demonstrate two workflows in a SARS-CoV-2 model system for simultaneous detection of ag-specific B and T cells within the same sample using **(1)** Dextramer® and Klickmer® reagents in combination with flow cytometry or **(2)** dCODE Dextramer® and dCODE Klickmer® reagents in combination with the 10x Single-Cell Analysis System.


Workflows for simultaneous investigation of antigen-specific B and T cells in blood samples



Dextramer® and Klickmer® reagents reveal changes in magnitude and kinetics of antigen-specific B and T cells upon vaccination



dCODE Klickmer® reagents identify BCR clonotypes



Conclusion

We have demonstrated two workflows for the simultaneous detection and characterization of ag-specific B and T cells in blood samples. **(1) Workflow 1** combines Dextramer® and Klickmer® reagents with flow cytometry to detect and characterize ag-specific B cells, CD4⁺ and CD8⁺ T cells. The workflow was demonstrated using SARS-CoV-2 as a model system to detect changes ag-specific B and T cells upon vaccination.

(2) Workflow 2 combines dCODE Dextramer® and dCODE Klickmer® with single-cell RNA seq to enable the examination of T and B cells in the same sample at the individual cell level and facilitates the evaluation of specific target BCR/TCR clonotypes. Although we are in the process of finalizing the sequencing of all samples, as an example we have presented data on BCRs for a single donor upon vaccination.