

# TCR Dextramer® reagents can measure antigen presentation on cells

Thomas Holberg Blicher, Kevin Lenogue, Liselotte Brix

Immudex ApS, Virum, Denmark


## Introduction

To successfully develop and apply TCR or T cell-based immunotherapies, the specificity and sensitivity of the selected TCR must first be validated before proceeding to clinical development. Here, the detection and quantification of antigen-presenting cells [APC] is important for **1)** stratification and selection of patients with demonstrated expression of the target antigen, **2)** confirming tissue-specific expression of the target antigen, and **3)** monitoring target expression during treatment. To support such efforts, we have developed high-avidity TCR Dextramer reagents to allow detection of peptide presentation by antigen-presenting cells. This study presents the use of such reagents as an analytical tool for evaluating target expression on the cell surface of antigen-presenting cells.

## Conclusions

- TCR Dextramer® reagents can be used in a simple workflow to detect peptide presentation at the surface of antigen-presenting cells, here shown for peptide-pulsed T2 cells and PBMC samples.
- Peptide presentation is detectable  $\geq 2$  nM peptide-pulsing concentration on T2 cells,  $\geq 20$  nM on PBMC samples.
- TCR Dextramer® detection sensitivity on T2 cells is only slightly influenced by the presence of irrelevant peptides.
- TCR Dextramer® is ideally suited to develop novel techniques for the detection of antigen presenting cells.

## TCR Dextramer reagents




TCR monomers are produced in *E. coli*, refolded and attached to a fluorescent Dextramer® backbone. Antigen-presenting cells [APC] can be stained with TCR Dextramer® reagents like conventional pMHC Dextramer® reagents on T cells.

## Experimental procedure


### TCRs used in this study

TCR	TCR Specificity	$K_d$
A	HLA-A*02:01/SLLMWITQV	48 pM
B	HLA-A*02:01/SLLMWITQV	32 pM



Peptide-pulsed cells are stained with TCR Dextramer® reagents based on high and low affinity TCRs. The mean fluorescence intensity of the relevant cell population [lymphocytes/monocytes] is measured using flow cytometry.

## Staining of peptide-pulsed cells using TCR Dextramer reagents



A) Peptide titration on T2 cells, B) TCR Dextramer® titration on peptide-pulsed T2 cells, C) Peptide titration on peptide-pulsed T2 cells using either target peptide alone or target peptide supplemented with a control peptide to a constant total peptide concentration of 2  $\mu$ M. Target:control peptide ratio is indicated on the graph. D) Peptide titration on PBMC sample; data shown for monocyte population. Target peptide: SLLMWITQV,  $K_d$  [HLA-A2] = 5.6 nM; control peptide: ALIAPVHAV,  $K_d$  [HLA-A2] = 7.0 nM.