

# A workflow for detection of antigen-specific HLA-E restricted T cells in PBMC from healthy individuals


Lyn Zhu, Dilek Inekci, Liselotte Brix  
Immudex ApS, Virum, Denmark




## Introduction

The non-classical human HLA-E molecule exhibits monomorphic functionality and presents a promising target for immunotherapeutic strategies against cancer and infectious diseases. A reliable assay for detecting HLA-E restricted antigen-specific T cells is essential for advancing research in this field. However, due to the inherent instability of HLA-E, both its production and the development of HLA-E specific detection reagents are challenging. Moreover, studies have demonstrated that HLA-E can be recognized by both NK cells and T cells via the CD94/NKG2 receptors, highlighting the importance of blocking these non-TCR-mediated interactions and implementing a specific gating strategy to accurately characterize antigen-specific HLA-E restricted CD8+ T cells by flow cytometry-based assays. Here, we present a workflow to detect such T cells in PBMCs from healthy individuals by using HLA-E Dextramer® technology and an optimized staining protocol.


## Workflow for detection of antigen-specific HLA-E restricted T cells



## Gating strategy



## Results



## Conclusions

- We were able to make stable HLA-E/peptide Dextramer® reagents using MHCp monomers made by refolding or a peptide receptive HLA-E monomer.
- By using anti-CD94 antibody we were able to prevent the binding of the HLA-E to the CD94/NKG2 receptors, allowing specific detection of TCR-mediated binding of T cells.
- Relevant cell subsets were identified by the designed antibody panel and gating strategy.
- Antigen-specific HLA-E restricted CD8+ T cells were identified by dual-staining of HLA-E Dextramer® reagents.
- we have established an optimized workflow for detection of HLA-E restricted antigen-specific T cells in PBMCs.



MM0278.01  
For research use only. Not for use in diagnostic or therapeutic procedures.  
© Copyright Immudex ApS, Denmark, 2024. All rights reserved.

IMMUDEX®