

# Developing a Platform for Screening T cell receptor Specificity and Ranking of Affinity

Amalie S. Schøller, Thomas H. Blicher and Liselotte Brix  
Immudex ApS, Virum, Denmark



DOWNLOAD NOW

FOLLOW US ON LinkedIn

## Introduction


Highly potent T cell therapy relies on the specificity as well as the affinity of the T cell receptor towards its MHC:peptide complex (MHCp). Here we present a platform that enables both production and functional validation of soluble TCR and TCR mimic (TCRm) candidates. By exploiting the high sensitivity of Dextramer® reagents, this flexible platform also allows medium to high-throughput screening for fine specificity and cross-reactivity. The advantages of this is furthermore that the binding characteristics for each TCR or TCRm Dextramer® reagent can be transformed into a binding score, which enables ranking of the candidates based on their affinity and overall avidity for MHCp complexes.

## Conclusions

We have established a platform for production and functional validation of TCR and TCRm molecules of choice, which also enables:


- Production of functional and sensitive TCR and TCRm Dextramer® reagents.
- Ranking of TCR molecules based on their affinity for any MHCp in an artificial cell system using flow-cytometry.
- A flexible platform for high-throughput screening for fine-specificity and cross-reactivity of TCR and TCRm candidates.

## Platform for Production and Validation




Separate chains of the monomer of choice, TCR or TCRm, were produced in *E. coli* and refolded *in vitro*. The resulting TCRm monomers were multimerized by attaching them to a fluorescent, PE-labelled Dextramer® backbone and finally tested for binding to MHCp-coated beads carrying either cognate or non-cognate MHCp complexes. The beads thus behave as an artificial cell display for validation of TCRm : MHCp binding.

## Experimental procedure



**A)** Heavy and light chains were successfully expressed in *E. coli* for three TCRm's. **B)** The expressed chains were purified and analysed by SDS PAGE demonstrating bands corresponding to the expected size of the light (L) and heavy (H) chains for all three TCRm constructs. **C)** The light and heavy chains were refolded and purified by FPLC; Fractions from the purification of the refolded TCRms were analyzed by SDS PAGE. **TCRm1** is shown as an example. Correct stoichiometry of heavy and light chain in the refolded **TCRm1** is shown by the equal intensity of the two bands corresponding to the light and heavy chains of **TCRm1** under reducing conditions. **D)** Successful recognition and binding of TCRm Dextramer® reagent to the target MHCp vs. binding to non-cognate MHCp complexes was measured as the increase in fluorescence intensity of PE.

## TCR affinity ranking



By varying the number of TCR molecules on the TCR Dextramer®, the TCR Dextramer molecules binds to the MHCp-loaded beads bind in a manner, which depends on the TCR:MHCp monomer affinity.

**A)** Binding of fluorescent TCR Dextramer reagents (T cell surrogate) to MHCp-loaded beads (target cell surrogate) was assessed using flow cytometry. Binding characteristics of each TCR Dextramer reagent was transformed into a binding score (R) by taking the ratio of the bead (geometric mean fluorescence intensity (gMFI) values at points A and B. **B)** The binding score correlates with TCR:MHCp equilibrium dissociation constants, and the empirical fit can be used to assess affinity directly or to simply perform affinity ranking. By choosing the points A and B, the affinity ranking can be adjusted for maximal sensitivity at either low, medium, or high affinity.

## Perspectives

This platform enables affinity assessment of TCR and TCRm in a high throughput format

- Assessment of fine specificity by screening against an array of MHCp complexes carrying variants of the target epitope
- Cross reactivity assessment by screening specific TCR or TCR candidates against an array of different/unrelated pMHC complexes
- Affinity ranking of chosen candidates by targeting the same or different MHCp complexes

