





| | MHC Dextramer® | | dCODE Dextramer® [HiT] - GOLD | | dCODE Dextramer® [10x compatible] - GOLD | | dCODE Dextramer® [RiO] - GOLD | | |
|----------------------------|--|--|---|-------------|--|--------------|--|--------------|--|
| |  | |  | |  | |  | | |
| Product description | Fluorochrome conjugated MHC Dextramer® | | DNA oligo barcoded MHC Dextramer® | | DNA oligo barcoded Dextramer® compatible with 10x Genomics Single-Cell Immune Profiling Solution | | DNA oligo barcoded Dextramer® compatible with BD Rhapsody™ Single-Cell Analysis System | | |
| Techniques used | Flow cytometry | | PCR and Next-Generation-Sequencing | | Multi-omics | | Multi-omics | | |
| Intended use | Identify virus-specific T cells in blood | | High-throughput, High-plex analysis of virus-specific T cells in blood | | Single-cell analysis of antigen-specific populations; gene expression; TCR sequencing | | Single-cell analysis of antigen-specific populations; gene expression | | |
| Label | APC, PE, FITC or none | | Up to 200 different DNA barcodes | | Up to 1000 different DNA barcodes | | Up to 1000 different DNA barcodes | | |
| Test Size / Reagent | 50 or 150 tests | | 25, 50 or 150 tests | | 25, 50 or 150 tests | | 25, 50 or 150 tests | | |
| SARS-CoV-2 Panel | XL | S | XL | S | XL | S | XL | S | |
| Catalog Number | Spike | COV1 XL | COV1 | dC01-COV1XL | dC01-COV1 | dCX01-COV1XL | dCX01-COV1 | dCR01-COV1XL | dCR01-COV1 |
| | A*0201 | COV2 XL | COV2 | dC01-COV2XL | dC01-COV2 | dCX01-COV2XL | dCX01-COV2 | dCR01-COV2XL | dCR01-COV2 |
| | Multi Allele | COV3 XL | COV3 | dC01-COV3XL | dC01-COV3 | dCX01-COV3XL | dCX01-COV3 | dCR01-COV3XL | dCR01-COV3 |
| | Custom | No Cat. No. Please specify number of Dextramer® reagents and which HLA-allele + epitope combination(s) you want. | | | | | | | |
| Price | Please contact us at customer@immudex.com with information about which Panel you are interested in, the label, the number of tests, the catalog number, and we will come back to you with information about the price and the estimated delivery time. | | | | | | | | |
| How to order | | | | | | | | | |
| Questions | | | | | | | | | We are happy to answer any questions you may have. |

Immudex is the sole manufacturer and provider of dCODE Dextramer® reagents, and support related to these products is through Immudex. For research use only. Not for use in diagnostic or therapeutic procedures.