





	MHC Dextramer®		dCODE Dextramer® [HiT] - GOLD		dCODE Dextramer® [10x compatible] - GOLD		dCODE Dextramer® [RiO] - GOLD		
									
Product description	Fluorochrome conjugated MHC Dextramer®		DNA oligo barcoded MHC Dextramer®		DNA oligo barcoded Dextramer® compatible with 10x Genomics Single-Cell Immune Profiling Solution		DNA oligo barcoded Dextramer® compatible with BD Rhapsody™ Single-Cell Analysis System		
Techniques used	Flow cytometry		PCR and Next-Generation-Sequencing		Multi-omics		Multi-omics		
Intended use	Identify virus-specific T cells in blood		High-throughput, High-plex analysis of virus-specific T cells in blood		Single-cell analysis of antigen-specific populations; gene expression; TCR sequencing		Single-cell analysis of antigen-specific populations; gene expression		
Label	APC, PE, FITC or none		Up to 200 different DNA barcodes		Up to 1000 different DNA barcodes		Up to 1000 different DNA barcodes		
Test Size / Reagent	50 or 150 tests		25, 50 or 150 tests		25, 50 or 150 tests		25, 50 or 150 tests		
SARS-CoV-2 Panel	XL	S	XL	S	XL	S	XL	S	
Catalog Number	Spike	RX07	RX06	RX13	RX12	RX19	RX18	RX25	RX24
	A*0201	RX09	RX08	RX15	RX14	RX21	RX20	RX27	RX26
	Multi Allele	RX11	RX10	RX17	RX16	RX23	RX22	RX29	RX28
	Custom	No Cat. No. Please specify number of Dextramer® reagents and which HLA-allele + epitope combination(s) you want.							
Price	Please contact us at customer@immudex.com with information about which Panel you are interested in, the label, the number of tests, the catalog number, and we will come back to you with information about the price and the estimated delivery time.								
How to order									
Questions									We are happy to answer any questions you may have.

Immudex is the sole manufacturer and provider of dCODE Dextramer® reagents, and support related to these products is through Immudex. For research use only. Not for use in diagnostic or therapeutic procedures.