

Accelerate Your Immuno-Oncology Research from Early Discovery to Clinical Trials

Version 1. October 2021

Recent advances in revealing how the immune system can recognize and target cancer cells have led to novel immunotherapeutic strategies within the immune-oncology field. Although immunotherapy has brought new and improved treatments to patients resistant to the classical ones, scientists worldwide continue their research to develop more potent personalized cancer vaccines with minimal toxicity effect.

IMPORTANCE OF DEFINING SUITABLE TUMOR ANTIGENS

When tumor cells are formed, they express antigens on the cell surface as part of the major histocompatibility complex (MHC) or release them into the bloodstream. The immune system can recognize these tumor antigens, and that triggers an immune response against cancer. Thus, specific tumor antigens can be valuable markers in cancer diagnostics and used as tumor antigen vaccines in personalized cancer immunotherapy.


Tumor antigens are generally classified into two categories based on their expression pattern:

- **tumor-specific antigens (TSA)** primarily found in cancer cells. TSA are mutation-derived neoantigens and can cause the tumor due to mutation in the conserved gene (p53 genes).
- **tumor-associated antigens (TAA)** highly expressed on tumor cells but can also be found at lower levels on healthy cells. TAA are produced due to the mutation in genes unrelated to tumor formation.

Other tumor antigens include tissue differentiation antigens, cancer-germline antigens, oncofetal antigens, etc. Thus, defining suitable tumor antigens is essential for triggering efficient and durable host immune responses against cancer.


MOVE FORWARD WITH YOUR CANCER RESEARCH STRATEGY

Immudex Dextramer® technology can help you detect, quantify or isolate cancer-specific CD8+ and CD4+ T cells in peripheral blood, tumor biopsies, or tissue without risking missing important information.



SOLUTIONS THAT ADAPT TO YOUR RESEARCH NEEDS

Immudex is an established expert in cancer-specific immune monitoring from early discovery to clinical trials. Our Dextramer® and dCODE Dextramer® technologies can help you characterize cancer immunity across different platforms, allowing you to continue directly from in-situ, to flow cytometry, and move onto NGS or single-cell multi-omics.



- **MHC I** and **MHC II Dextramer®** reagents comprise a broad selection of T-cell epitopes to monitor the whole spectrum of cancer-specific T-cell responses by flow cytometry, including T cells with rare and low-affinity receptors. Hand-picked and **ready-to-use melanoma** and **cancer-testis panels** are also available for your research.
- **MHC I Dextramer® In-Situ Staining** reagents can help detect distinct localization, distribution, and abundance of cancer-specific CD8+ T cells within diverse tissue niches in health and disease.
- **dCODE Dextramer®** reagents for epitope discovery or T-cell profiling in bulk solution. The unique DNA barcode enables the detection of antigen-specific T cells by PCR followed by next generation sequencing (NGS) and single-cell multi-omics.
- **U-Load Dextramer®** reagents to easily make your choice of customized multimers directly in the lab from a broad range of peptide-receptive MHC I and MHC II alleles and your own peptide, via flow cytometry or single-cell multi-omics.
- **Klickmer®** reagents can be customized to help you to explore immunity beyond T cells by looking into other immune cells like B cells or study TCR-ligand interactions in cancer cells by flow cytometry or single-cell multi-omics.

Learn more or contact us to discuss your research need at customer@immudex.com

SELECT THE MOST SUITABLE TUMOR-SPECIFIC ANTIGEN FOR YOUR RESEARCH

To support clinicians and scientists worldwide in their cancer research, we want to share a comprehensive list of human cancer antigens recognized by T cells. This list of cancer-specific T-cell epitopes reflects the research of many dedicated scientists published in peer-reviewed papers. You can now select the epitopes from the list that best fit your cancer research.

Published Human Cancer T-Cell Epitopes

Version 1. October 2021

HLA Allele	Peptide	Antigen	Disease	Reference
A*0101	TLDTLTAFY	MSLN	Cancer	3
A*0101	FTELTGFEF	Survivin	Cancer	4
A*0101	TSEKRPFMCAY	WT1	Cancer	7, 8
A*0101	EADPTGHSY	MAGE-1	Melanoma	2
A*0101	EVDPIGHLY	MAGE-3	Melanoma	1, 2, 10
A*0101	KSDICTDEY	Tyrosinase	Melanoma	9
A*0101	SSDYVIPIGTY	Tyrosinase	Melanoma	11
A*0201	QQAHLWCV	ABL1	Cancer	13
A*0201	ALDVYNGLL	ACPP	Cancer	12
A*0201	GLQHWVPEL	BA46	Cancer	17
A*0201	NLFETPVEA	BA46	Cancer	17
A*0201	KLDVGNAEV	BAP31	Cancer	14
A*0201	ALSPVPPVV	BCL-2	Cancer	15
A*0201	WLSLKTLLSL	BCL-2	Cancer	16, 24
A*0201	YLNRLHTWI	BCL-2	Cancer	18
A*0201	YLNDHLEPWI	BCL-X	Cancer	16, 19
A*0201	GVRGRVEEI	BCR-ABL	Cancer	20
A*0201	CLPSPSTPV	BMI1	Cancer	21, 22
A*0201	TLQDIVYKL	BMI1	Cancer	21, 22
A*0201	MLMAQEALAFI	CAMEL	Cancer	23
A*0201	ALYLMELTM	CB9L2	Cancer	24, 25
A*0201	YLISGDSPV	CD33	Cancer	26
A*0201	YLSGANLNL	CEACAM	Cancer	27, 28
A*0201	ILGVLTSLV	DLK1	Cancer	29
A*0201	TLADFDPRV	EphA2	Cancer	30, 31
A*0201	SQADALKYV	EZH2	Cancer	32, 33

Published Human Cancer T-cell Epitopes

Version 1. October 2021

HLA Allele	Peptide	Antigen	Disease	Reference
A*0201	YMCSFLFNL	EZH2	Cancer	32
A*0201	TLFWLLLT	VEGFR1	Cancer	149, 152
A*0201	EIWTHSYKV	FOLR1	Cancer	34
A*0201	FVGEFFTDV	Glypican 3	Cancer	35, 36, 145, 146
A*0201	YLEPGPVT	gp100	Cancer	37, 38
A*0201	QLFEELQEL	HO-1	Cancer	39
A*0201	LLGLPLGL	Heparanase	Cancer	40, 41
A*0201	ILHDGAYSL	HER2	Cancer	42
A*0201	KIFGSLAFL	HER2	Cancer	43, 44
A*0201	LIAHNQVRQV	HER2	Cancer	45
A*0201	ILSLELMKL	HMMR	Cancer	46, 47
A*0201	ALPFGFILV	IL13Ra	Cancer	30, 153
A*0201	ALMEQQHYV	ITGB8	Cancer	48
A*0201	VISNDVCAQV	KLK	Cancer	49
A*0201	FIYDFCIFGV	Lengsin	Cancer	50, 51
A*0201	QLCPICRAPV	LIVIN	Cancer	52, 53
A*0201	LLASIAAGL	LY6K	Cancer	149
A*0201	SLLFLLFSL	MSLN	Cancer	54, 55, 173
A*0201	VLPLTVAEV	MSLN	Cancer	54, 55, 173
A*0201	AQCQETIRV	Midkine	Cancer	56
A*0201	SLFLGILSV	MS4A1	Cancer	57
A*0201	SLLMWITQC	NY-ESO-1	Cancer	58, 59, 60
A*0201	SLLMWITQV	NY-ESO-1	Cancer	60, 61
A*0201	YLYQWLGA	BGLAP	Cancer	154
A*0201	GLAPPQHLIRV	p53	Cancer	62, 63

Published Human Cancer T-Cell Epitopes

Version 1. October 2021

HLA Allele	Peptide	Antigen	Disease	Reference
A*0201	KLCPVQLWV	p53	Cancer	62, 63
A*0201	LLGRNSFEV	p53	Cancer	62, 63
A*0201	RMPEAAPPV	p53	Cancer	62, 63
A*0201	SLPPPGRV	p53	Cancer	62, 63
A*0201	STPPPGRV	p53	Cancer	64
A*0201	YLGSYGFRL	p53	Cancer	62, 63
A*0201	ELSDSLGPV	PASD1	Cancer	65
A*0201	QLLDGFMITL	PASD1	Cancer	66
A*0201	SIDWFMVTV	PLAC1	Cancer	67, 68
A*0201	VLQELNVTV	Pr1	Cancer	69, 70
A*0201	ALYVDSLFFL	PRAME	Cancer	71, 72, 73
A*0201	VLDGLDVLL	PRAME	Cancer	71, 72, 73
A*0201	YLQWIEFSI	Prominin1	Cancer	74
A*0201	FLTPKKLQCV	PSA	Cancer	75
A*0201	KLQCVDLHV	PSA	Cancer	76
A*0201	AILALLPAL	PSCA	Cancer	77, 78
A*0201	MMNDQLMFL	PSMA	Cancer	79, 80
A*0201	VLAGGFLL	PSMA	Cancer	79, 80
A*0201	ALWPWLLMAT	RNF43	Cancer	81
A*0201	RLAEYQAYI	SART3	Cancer	82, 83
A*0201	MLAVFLPIV	STEAP1	Cancer	84
A*0201	TLPPAWQPFL	Survivin	Cancer	5, 85
A*0201	LMLGEFLKL	Survivin	Cancer	24, 86
A*0201	LTLGEFLKL	Survivin-3a	Cancer	87

Published Human Cancer T-cell Epitopes

Version 1. October 2021

HLA Allele	Peptide	Antigen	Disease	Reference
A*0201	YLIELIDRV	TACE	Cancer	88
A*0201	FLFLRNFSL	TARP	Cancer	89, 90
A*0201	FLPSPLFFFL	TARP 2M	Cancer	89, 90
A*0201	ALLTSRLRFI	Telomerase	Cancer	91
A*0201	ILAKFLHWL	Telomerase	Cancer	92, 93, 94
A*0201	YLQVNSLQTV	Telomerase	Cancer	95, 96
A*0201	RLSSCVPVA	TGFβ	Cancer	97, 98
A*0201	FLYDDNQRV	topII	Cancer	14
A*0201	LLLLTVLTV	Mucin	Cancer	99, 100
A*0201	RMFPNAPYL	WT1	Cancer	101, 102, 103, 106
A*0201	SLGEQQYSV	WT1	Cancer	103, 104, 106
A*0201	VLDFAAPGA	WT1	Cancer	105, 106
A*0201	IMDQVPFSV	gp100	Melanoma	37, 107
A*0201	ITDQVPFSV	gp100	Melanoma	108, 109
A*0201	KTWGQYWQV	gp100	Melanoma	108, 109, 148
A*0201	GLYDGMEHL	MAGE-10	Melanoma	110, 111
A*0201	FLWGPRALV	MAGE-3	Melanoma	112, 113
A*0201	GVYDGREHTV	MAGE-4	Melanoma	114, 115, 116
A*0201	KVLEYVIKV	MAGE-A1	Melanoma	117, 118, 119
A*0201	YLQLVFGIEV	MAGE-A2	Melanoma	120, 121
A*0201	KVAELVHFL	MAGE-A3	Melanoma	111, 122, 123, 124
A*0201	FLAMLKNTV	MAGE-C1	Melanoma	125, 126
A*0201	ELAGIGILTV	MART-1	Melanoma	120, 127, 128
A*0201	YMDGTMSQV	Tyrosinase	Melanoma	37, 129, 131

Published Human Cancer T-cell Epitopes

Version 1. October 2021

HLA Allele	Peptide	Antigen	Disease	Reference
A*0201	LLVPTCVFLV	Endosialin	Cancer	29, 191
A*0201	GMLGMVSGL	NRP-1	Cancer	29, 191
A*0201	SLFEPPPPG	PSMA	Cancer	192
A*0301	RIAAWMATY	BCL-2L1	Cancer	24
A*0301	ALLAVGATK	gp100	Cancer	132, 133
A*0301	RLGLQVRKNK	RhoC	Cancer	24, 134
A*0301	RISTFKNWPK	Survivin-3a	Cancer	4, 135, 136
A*0301	YMVPFIPLYR	Tyrosinase	Melanoma	137
A*2402	DYLQYVLQI	BCL-2A1	Cancer	138, 139
A*2402	EYRALQLHL	Carbonic anhydrase	Cancer	140
A*2402	EYYELFVNI	DEP DC1	Cancer	141, 142, 143
A*2402	IYTWIEDHF	FOX M1	Cancer	142, 144, 149, 150
A*2402	EYILSLEEL	Glypican 3	Cancer	145, 146
A*2402	VYFFLPDHL	gp100	Cancer	147, 148
A*2402	KWLISPVKI	HJURP	Cancer	149, 150
A*2402	VYLRVRPLL	KIF20A	Cancer	151, 160, 161
A*2402	RYCNLEGPII	LY6K	Cancer	155, 156
A*2402	EYCPGGNLF	MELK	Cancer	149, 150
A*2402	VYGIRLEHF	CDCA1	Cancer	157, 158, 159
A*2402	DYLNEWGSRF	p-Cadherin	Cancer	151, 160, 161
A*2402	CYASGWGSI	PSA	Cancer	162, 163
A*2402	AYACNTSTL	Survivin	Cancer	87, 164, 165, 166
A*2402	SYRNEIAYL	TTK	Cancer	167, 168
A*2402	CYTWNQMNL	WT1	Cancer	169, 170, 171

Published Human Cancer T-cell Epitopes

Version 1. October 2021

HLA Allele	Peptide	Antigen	Disease	Reference
A*2402	AFLPWHRLF	Tyrosinase	Melanoma	172, 173
A*6801	TVSGNILTIR	NY-ESO-1	Cancer	174
B*0702	APRGVRMAV	CAMEL	Cancer	175
B*0702	TPNQRQNVC	P2X5a	Cancer	176, 177
B*0702	RVRFFFPSL	MAGE-A1	Melanoma	178, 179
B*0702	LPWHRLFLL	Tyrosinase	Melanoma	137
B*3501	MPFATPMEA	NY-ESO-1	Cancer	180, 181, 182, 184
B*3501	EAAGIGILTY	MART-1	Melanoma	183
B*3501	MPFATPMEAEL	NY-ESO-1	Melanoma	10, 184
B*5101	MPFATPMEA	NY-ESO-1	Cancer	185
DPB1*0401	SLLMWITQCFLPVF	NY-ESO	Cancer	186, 187, 188
DPB1*0401	TQHFVQENYLEY	MAGE-A3	Cancer	189, 190

REFERENCES

1. Raman, M., Rizkallah, P., Simmons, R., *et al.* (2016). Direct molecular mimicry enables off-target cardiovascular toxicity by an enhanced affinity TCR designed for cancer immunotherapy. <https://doi.org/10.1038/srep18851>
2. Karanikas, V., *et al.* (2003). Monoclonal Anti-MAGE-3 CTL Responses in Melanoma Patients Displaying Tumor Regression after Vaccination with a Recombinant Canarypox Virus. <https://doi.org/10.4049/jimmunol.171.9.4898>
3. Laheru, D., *et al.* (2008). Allogeneic Granulocyte Macrophage Colony-Stimulating Factor-Secreting Tumor Immunotherapy Alone or in Sequence with Cyclophosphamide for Metastatic Pancreatic Cancer: A Pilot Study of Safety, Feasibility, and Immune Activation. doi: [10.1158/1078-0432.CCR-07-0371](https://doi.org/10.1158/1078-0432.CCR-07-0371)
4. Reker, S., *et al.* (2004). Identification of novel survivin-derived CTL epitopes. doi: [10.4161/cbt.3.2.611](https://doi.org/10.4161/cbt.3.2.611)
5. Schmitz, M., *et al.* (2000). Generation of survivin-specific CD8 β T effector cells by dendritic cells pulsed with protein or selected peptides.
6. Reker, S., *et al.* (2004). HLA-B35-restricted immune responses against survivin in cancer patients. doi: [10.1002/ijc.11634](https://doi.org/10.1002/ijc.11634)
7. Asemissen, A. M., *et al.* (2006). Identification of a Highly Immunogenic HLA-A*01-Binding T Cell Epitope of WT1. doi: [10.1158/1078-0432.CCR-06-1337](https://doi.org/10.1158/1078-0432.CCR-06-1337)

REFERENCES

8. Weber, G., *et al.* (2009). WT1 peptide-specific T cells generated from peripheral blood of healthy donors: possible implications for adoptive immunotherapy after allogeneic stem cell transplantation. <https://doi.org/10.1038/leu.2009.70>
9. Kittlesen, D. J., *et al.* (1998). Human Melanoma Patients Recognize an HLA-A1-Restricted CTL Epitope from Tyrosinase Containing Two Cysteine Residues: Implications for Tumor Vaccine Development.
10. Sahin, U., *et al.* (2020). An RNA vaccine drives immunity in checkpoint-inhibitor-treated melanoma. doi: <https://doi.org/10.1038/s41586-020-2537-9>
11. Kawakami, Y., *et al.* (1998). Identification of new melanoma epitopes on melanosomal proteins recognized by tumor infiltrating T lymphocytes restricted by HLA-A1, -A2, and -A3 alleles. *J Immunol.*
12. Olson, B. M., *et al.*, (2010). HLA-A2-restricted T-cell epitopes specific for prostatic acid phosphatase. *Cancer Immunol Immunother.* 2010 Jun;59(6):943-53. doi: [10.1007/s00262-010-0820-6](https://doi.org/10.1007/s00262-010-0820-6)
13. Chiarella, P., *et al.* (2012). BCR/ABL1 Fusion Transcripts Generated from Alternative Splicing: Implications for Future Targeted Therapies in Ph+ Leukaemias. doi: <https://doi.org/10.2174/156652412800619996>
14. Penwell, A., *et al.* (2012). Development and validation of an HPLC/UV assay for separation and quantification of peptide antigens from a liposomal vaccine delivery platform. doi: <https://doi.org/10.1016/j.jpba.2012.03.046>
15. Wang, M., *et al.* (2007). Identification of an HLA-A*0201 restricted Bcl2-derived epitope expressed on tumors.
16. Jørgensen, N. G., *et al.* (2016). Peptide vaccination against multiple myeloma using peptides derived from anti-apoptotic proteins: a phase I trial. doi: [10.21037/sci.2016.11.09](https://doi.org/10.21037/sci.2016.11.09)
17. Jurtz, V. I., *et al.* (2018). NetTCR: sequence-based prediction of TCR binding to peptide-MHC complexes using convolutional neural networks.
18. Andersen, M. H., *et al.* (2005). Immunogenicity of Bcl-2 in patients with cancer. doi: <https://doi.org/10.1182/blood-2004-07-2548>
19. Sørensen, R. B., *et al.* (2007). Efficient tumor cell lysis mediated by a Bcl-X(L) specific T cell clone isolated from a breast cancer patient. doi: <https://link.springer.com/article/10.1007%2Fs00262-006-0200-4#citeas>
20. Volpe, G., *et al.* (2007). Alternative BCR/ABL Splice Variants in Philadelphia Chromosome-Positive Leukemias Result in Novel Tumor-Specific Fusion Proteins that May Represent Potential Targets for Immunotherapy Approaches.
21. Hughes, A., *et al.* (2017). CML patients with deep molecular responses to TKI have restored immune effectors and decreased PD-1 and immune suppressors. doi: <https://doi.org/10.1182/blood-2016-10-745992>
22. Agnes Yong, S. M., *et al.* (2009). In Chronic Myeloid Leukemia Cytotoxic T Cell Responses to BMI-1 Protein Correlate with Higher Expression of *BMI-1* in Leukemia Progenitors, and May Contribute to Improved Outcome After Allogeneic Stem Cell Transplantation. doi: <https://doi.org/10.1182/blood.V114.22.192.192>
23. Aarnoudse, C. A., *et al.* (1999). Interleukin-2-induced, melanomaspecific T cells recognize CAMEL, an unexpected translation product of LAGE-1. doi: [https://doi.org/10.1002/\(SICI\)1097-0215\(19990730\)82:3<442::AID-IJC19>3.0.CO;2-Z](https://doi.org/10.1002/(SICI)1097-0215(19990730)82:3<442::AID-IJC19>3.0.CO;2-Z)
24. Junker, N., *et al.* (2012). Tumor associated antigen specific T-cell populations identified in ex vivo expanded TIL cultures.

REFERENCES

25. Sørensen, R. B., *et al.* (2009). CD8 T-cell Responses against Cyclin B1 in Breast Cancer Patients with Tumors Overexpressing p53.
26. Bae, J., *et al.* (2004). Heteroclitic CD33 Peptide With Enhanced Anti-Acute Myeloid Leukemic Immunogenicity. doi: [10.1158/1078-0432.CCR-04-0322](https://doi.org/10.1158/1078-0432.CCR-04-0322)
27. Fauquembergue, E., *et al.* (2010). HLA-A*0201-restricted CEA-derived Peptide CAP1 Is Not a Suitable Target for T-cell-based Immunotherapy.
28. Mellstedt, H., *et al.* (2017). DNA immunization targeting carcinoembryonic antigen in colorectal cancer patients.
29. Zhao, X., *et al.* (2011). Intratumoral IL-12 Gene Therapy Results in the Crosspriming of Tc1 Cells Reactive Against Tumor-associated Stromal Antigens. doi: <https://doi.org/10.1038/mt.2010.295>
30. Pollack, I. F., *et al.* (2016). Immune responses and outcome after vaccination with glioma-associated antigen peptides and poly-ICLC in a pilot study for pediatric recurrent low-grade gliomas. doi: <https://doi.org/10.1093/neuonc/now026>
31. Hatano, M., *et al.* (2005). EphA2 as a Glioma-Associated Antigen: A Novel Target for Glioma Vaccines. doi: <https://doi.org/10.1593/neo.05277>
32. Steele, J., *et al.* (2006). The polycomb group proteins, BMI-1 and EZH2, are tumour-associated antigens. doi: <https://www.nature.com/articles/6603369#citeas>
33. Itoh, Y., *et al.* (2007). New peptides of the polycomb group protein enhancer of zeste homolog 2 with the potential to induce cancer-reactive cytotoxic T lymphocytes in human leukocyte antigen-A2+ prostate cancer patients.
34. Brown, T. A., *et al.* (2019). Final analysis of a phase I/IIa trial of the folate-binding protein-derived E39 peptide vaccine to prevent recurrence in ovarian and endometrial cancer patients. doi: <https://doi.org/10.1002/cam4.2378>
35. Yoshikawa T., *et al.* (2012). Abstract 3508: Zoledronate induce expansion of glypican-3 (GPC3) peptide specific cytotoxic T lymphocytes sufficient for adoptive cancer immunotherapy.
36. Yoshikawa T., *et al.* (2011). Abstract 2688: Establishment of high avidity CTL clones induced by HLA-A2 restricted glypican-3 (GPC3) peptide vaccine.
37. Yi, X., *et al.* (2018). Identification of the Risk HLA-A Alleles and Autoantigen in Han Chinese Vitiligo Patients and the Association of CD8+T Cell Reactivity with Disease Characteristics. doi: [10.12659/MSM.910515](https://doi.org/10.12659/MSM.910515)
38. Tan, M. P., *et al.* (2017). Human leucocyte antigen class I-redirected anti-tumour CD4⁺ T cells require a higher T cell receptor binding affinity for optimal activity than CD8⁺ T cells. Clin Exp Immunol, 187: 124-137. doi: <https://doi.org/10.1111/cei.12828>
39. Andersen, M. H., *et al.* (2009). Identification of heme oxygenase-1-specific regulatory CD8⁺ T cells in cancer patients. doi: [10.1172/JCI38739](https://doi.org/10.1172/JCI38739)
40. Domschke, C., *et al.* (2010). Effects of distant metastasis and peripheral CA 15-3 on the induction of spontaneous T cell responses in breast cancer patients.
41. Sommerfeldt, N., *et al.* (2006). Heparanase: A New Metastasis-Associated Antigen Recognized in Breast Cancer Patients by Spontaneously Induced Memory T Lymphocytes.
42. Anticoli, S., *et al.* (2018). Engineered exosomes emerging from muscle cells break immune tolerance to HER2 in transgenic mice and induce antigen-specific CTLs upon challenge by human dendritic cells. doi: <https://link.springer.com/article/10.1007/s00109-017-1617-2>

REFERENCES

43. Krokhotin, A., *et al.* (2019). Computationally Guided Design of Single-Chain Variable Fragment Improves Specificity of Chimeric Antigen Receptors. doi: <https://doi.org/10.1016/j.omto.2019.08.008>
44. Kuznetsova, M., *et al.* (2017). Generation of populations of antigen-specific cytotoxic T cells using DCs transfected with DNA construct encoding HER2/neu tumor antigens epitopes.
45. Angelos, D., *et al.* (2015). Identification of a Novel Immunogenic HLA-A*0201-Binding Epitope of HER-2/neu with Potent Antitumor Properties. doi: [10.4049/jimmunol.181.1.146](https://doi.org/10.4049/jimmunol.181.1.146)
46. Willemen, Y., *et al.* (2016). The tumor-associated antigen RHAMM (HMMR/CD168) is expressed by monocyte-derived dendritic cells and presented to T cells. doi: [10.18632/oncotarget.12170](https://doi.org/10.18632/oncotarget.12170)
47. Greiner, J., *et al.* (2006). Expression of tumor-associated antigens in acute myeloid leukemia: implications for specific immunotherapeutic approaches. doi: <https://doi.org/10.1182/blood-2006-01-023127>
48. Shraibman, B., *et al.* (2019). Identification of Tumor Antigens Among the HLA Peptidomes of Glioblastoma Tumors and Plasma.
49. Wilkinson, R., *et al.* (2012). Human kallikrein 4 signal peptide induces cytotoxic T cell responses in healthy donors and prostate cancer patients. doi: <https://doi.org/10.1007/s00262-011-1095-2>
50. Nakatsugawa, M., *et al.* (2011). Identification of an HLA-A*0201-restricted cytotoxic T lymphocyte epitope from the lung carcinoma antigen, Lentsin. doi: <https://doi.org/10.3892/ijo.2011.1089>
51. Nakatsugawa, M., *et al.* (2009). Novel spliced form of a lens protein as a novel lung cancer antigen, Lentsin splicing variant 4. doi: <https://doi.org/10.1111/j.1349-7006.2009.01187.x>
52. Andersen, M. H., *et al.* (2004). The Melanoma Inhibitor of Apoptosis Protein: A Target for Spontaneous Cytotoxic T Cell Responses. doi: <https://doi.org/10.1046/j.0022-202X.2004.22242.x>
53. Sørensen, M. F., *et al.* (2009). Melanoma Inhibitor of Apoptosis Protein (ML-IAP) Specific Cytotoxic T Lymphocytes Cross-React with an Epitope from the Auto-Antigen SS56. doi: <https://doi.org/10.1038/jid.2009.10>
54. Klaver, Y., *et al.* (2019). Autologous dendritic cell therapy in mesothelioma patients enhances frequencies of peripheral CD4 T cells expressing HLA-DR, PD-1 or ICOS.
55. Miyazawa, M., *et al.* (2011). Dendritic cells adenovirally-transduced with full-length mesothelin cDNA elicit mesothelin-specific cytotoxicity against pancreatic cancer cell lines in vitro. doi: <https://doi.org/10.1016/j.canlet.2011.02.013>
56. Kerzerho, M., *et al.* (2010). The Angiogenic Growth Factor and Biomarker Midkine Is a Tumor-Shared Antigen. doi: <https://doi.org/10.4049/jimmunol.0901014>
57. Milcent, B., *et al.* (2019). Presence of T cells directed against CD20-derived peptides in healthy individuals and lymphoma patients. doi: <https://doi.org/10.1007/s00262-019-02389-7>
58. D'Angelo, S. P., *et al.* (2018). Antitumor Activity Associated with Prolonged Persistence of Adoptively Transferred NY-ESO-1 c259T Cells in Synovial Sarcoma. doi: [10.1158/2159-8290.CD-17-1417](https://doi.org/10.1158/2159-8290.CD-17-1417)
59. Ramachandran, I., *et al.* (2019). Systemic and local immunity following adoptive transfer of NY-ESO-1 SPEAR T cells in synovial sarcoma. doi: <https://doi.org/10.1186/s40425-019-0762-2>

REFERENCES

60. Shenderov, E., *et al.* (2021). Generation and characterization of HLA-A2 transgenic mice expressing the human TCR 1G4 specific for the HLA-A2 restricted NY-ESO-1₁₅₇₋₁₆₅ tumor-specific peptide. doi: [10.1136/jitc-2021-002544](https://doi.org/10.1136/jitc-2021-002544)
61. Neek, M., *et al.* (2018). Co-delivery of human cancer-testis antigens with adjuvant in protein nanoparticles induces higher cell-mediated immune responses. doi: <https://doi.org/10.1016/j.biomaterials.2017.11.022>
62. Svane, I. M., *et al.* (2006). Dendritic cells loaded with wild-type p53 peptides induce T cell immunity and disease stabilization in patients with advanced breast cancer.
63. Svane, I. M., *et al.* (2008). Alterations in p53-specific T cells and other lymphocyte subsets in breast cancer patients during vaccination with p53-peptide loaded dendritic cells and low-dose interleukin-2. doi: <https://doi.org/10.1016/j.vaccine.2008.06.085>
64. Würtzen, P. A., *et al.* (2002). A HLA-A2 restricted human CTL line recognizes a novel tumor cell expressed p53 epitope. *Int. J. Cancer*, 99: 568-572. doi: <https://doi.org/10.1002/ijc.10375>
65. Ait-Tahar, K., *et al.* (2009). Cytolytic T-cell response to the PASD1 cancer testis antigen in patients with diffuse large B-cell lymphoma. *British Journal of Haematology*, 146: 396-407. doi: <https://doi.org/10.1111/j.1365-2141.2009.07761.x>
66. Joseph-Pietras, D., *et al.* (2010). DNA vaccines to target the cancer testis antigen PASD1 in human multiple myeloma. doi: <https://doi.org/10.1038/leu.2010.196>
67. Taheri-Anganeh, M., *et al.* (2020). In silico Evaluation of PLAC1-fluC As a Chimeric Vaccine against Breast Cancer. doi: [10.29252/ibj.24.3.173](https://doi.org/10.29252/ibj.24.3.173)
68. Liu, W., *et al.* (2011). Identification of a novel HLA-A2-restricted cytotoxic T lymphocyte epitope from cancer-testis antigen PLAC1 in breast cancer.
69. Yong, A. S. M., *et al.* (2006). Cytotoxic T Lymphocyte Responses to PR1 Peptide in Chronic Myeloid Leukemia Patients Inversely Correlate with Proteinase 3 and Elastase Expression but Donor PR1 Responses Determine Survival after Stem Cell Transplantation. doi: <https://doi.org/10.1182/blood.V108.11.3697.3697>
70. Alatrash, H., *et al.* (2012). Multiple Myeloma Cross-Presents Neutrophil Elastase and Proteinase 3 and Becomes Susceptible to PR1-Targeting Immunotherapy. doi: <https://doi.org/10.1182/blood.V120.21.4973.4973>
71. Gejman, R. S., *et al.* (2018). Prospective identification of cross-reactive human peptide-MHC ligands for T cell receptor based therapies. doi: [10.1158/2326-6066.CIR-19-0745](https://doi.org/10.1158/2326-6066.CIR-19-0745)
72. Rezvani, K., *et al.* (2009). Ex vivo characterization of polyclonal memory CD8⁺ T-cell responses to PRAME-specific peptides in patients with acute lymphoblastic leukemia and acute and chronic myeloid leukemia. doi: <https://doi.org/10.1182/blood-2008-03-144071>
73. Yao Y, *et al.* (2013). Increased PRAME-Specific CTL Killing of Acute Myeloid Leukemia Cells by Either a Novel Histone Deacetylase Inhibitor Chidamide Alone or Combined Treatment with Decitabine. doi: <https://doi.org/10.1371/journal.pone.0070522>
74. Judkowski, J. J., *et al.* (2014). Identification of Novel Human Leukocyte Antigen-A*0201-Restricted, Cytotoxic T Lymphocyte Epitopes on CD133 for Cancer Stem Cell Immunotherapy. doi: <https://doi.org/10.5966/sctm.2013-0135>
75. Correale, P., *et al.* (1997). In vitro generation of human cytotoxic T lymphocytes specific for peptides derived from prostate-specific antigen.
76. Perambakam, S., *et al.* (2002). Induction of Tc2 cells with specificity for prostate-specific antigen from patients with hormone-refractory prostate cancer. doi: [10.1007/s00262-002-0281-7](https://doi.org/10.1007/s00262-002-0281-7)

REFERENCES

77. Feyerabend, S., *et al.* (2009). Novel multi-peptide vaccination in Hla-A2+ hormone sensitive patients with biochemical relapse of prostate cancer. doi: <https://doi.org/10.1002/pros.20941>
78. Kiessling, A., *et al.* (2002). Prostate stem cell antigen: Identification of immunogenic peptides and assessment of reactive CD8⁺ T cells in prostate cancer patients. doi: <https://doi.org/10.1002/ijc.10713>
79. Garetto, S., *et al.* (2009) Binding of prostate-specific membrane antigen to dendritic cells: a critical step in vaccine preparation. doi: [10.3109/14653240903164971](https://doi.org/10.3109/14653240903164971)
80. Vittes, G. E., *et al.* (2011). DNA fusion gene vaccines induce cytotoxic T-cell attack on naturally processed peptides of human prostate-specific membrane antigen. doi: <https://doi.org/10.1002/eji.201141518>
80. Vittes, G. E., *et al.* (2011). DNA fusion gene vaccines induce cytotoxic T-cell attack on naturally processed peptides of human prostate-specific membrane antigen. doi: <https://doi.org/10.1002/eji.201141518>
81. Hijikata, Y., *et al.* (2018). A phase I clinical trial of RNF43 peptide-related immune cell therapy combined with low-dose cyclophosphamide in patients with advanced solid tumors. doi: <https://doi.org/10.1371/journal.pone.0187878>
82. Noguchi, M., *et al.* (2018). Survival analysis of multiple peptide vaccination for the selection of correlated peptides in urological cancers. doi: <https://doi.org/10.1111/cas.13709>
83. Yoshiyama, K., *et al.* (2012). Personalized peptide vaccination in patients with refractory non-small cell lung cancer.
84. Stei, R. G., *et al.* (2019). Cognate Nonlytic Interactions between CD8⁺ T Cells and Breast Cancer Cells Induce Cancer Stem Cell-like Properties. doi: [10.1158/0008-5472.CAN-18-0387](https://doi.org/10.1158/0008-5472.CAN-18-0387)
85. Chang, D. Z., *et al.* (2007). Presence of survivin-specific cytotoxic T-lymphocytes (CTLs) in pancreatic cancer patients, and specific killing of human pancreatic carcinoma cells *in vitro* by survivin-specific CTLs.
86. Zhang, S., *et al.* (2020). Efficient induction of cytotoxic T lymphocytes in hepatocellular carcinoma using the HLA-A2-restricted survivin peptide *in vitro*. doi: <https://doi.org/10.1016/j.yexcr.2019.111741>
87. Karanikas, V., *et al.* (2008). Baseline levels of CD8⁺ T cells against survivin and survivin-2B in the blood of lung cancer patients and cancer-free individuals. doi: <https://doi.org/10.1016/j.clim.2008.07.024>
88. Sinnathamby, G., *et al.* (2011). ADAM metallopeptidase domain 17 (ADAM17) is naturally processed through major histocompatibility complex (MHC) class I molecules and is a potential immunotherapeutic target in breast, ovarian and prostate cancers. doi: <https://doi.org/10.1111/j.1365-2249.2010.04298.x>
89. Carlsson, B., *et al.* (2004). Generation of cytotoxic T lymphocytes specific for the prostate and breast tissue antigen TARP. doi: [10.1002/pros.20091](https://doi.org/10.1002/pros.20091)
90. Forsberg, O., *et al.* (2009). High frequency of prostate antigen-directed T cells in cancer patients compared to healthy age-matched individuals. doi: <https://doi.org/10.1002/pros.20858>
91. Inderberg-Suso, E-M., *et al.* (2012). Widespread CD4⁺ T-cell reactivity to novel hTERT epitopes following vaccination of cancer patients with a single hTERT peptide GV1001. doi: <https://doi.org/10.4161/onci.20426>

REFERENCES

92. Cole, D. K., *et al.* (2008). T cell receptor engagement of peptide-major histocompatibility complex class I does not modify CD8 binding. doi: <https://doi.org/10.1016/j.molimm.2007.12.009>
93. Huang, G., *et al.* (2012). Identification of a New Cytotoxic T-Cell Epitope p675 of Human Telomerase Reverse Transcriptase. doi: <https://doi.org/10.1089/cbr.2012.1193>
94. Cole, D. K., *et al.* (2017). Structural Mechanism Underpinning Cross-reactivity of a CD8⁺ T-cell Clone That Recognizes a Peptide Derived from Human Telomerase Reverse Transcriptase.
95. Brimnes, M.K., *et al.* (2012). Generation of autologous tumor-specific T cells for adoptive transfer based on vaccination, in vitro restimulation and CD3/CD28 dynabead-induced T cell expansion. doi: <https://doi.org/10.1007/s00262-011-1199-8>
96. Chu, C.S., *et al.* (2012). Phase I/II randomized trial of dendritic cell vaccination with or without cyclophosphamide for consolidation therapy of advanced ovarian cancer in first or second remission. doi: <https://doi.org/10.1007/s00262-011-1081-8>
97. Inderberg, E. M., *et al.* (2017). T cell therapy targeting a public neoantigen in microsatellite instable colon cancer reduces in vivo tumor growth.
98. Westdorp, H. Vaccine-induced multifunctional T cells recognize tumor-specific neoantigens in Lynch syndrome patients. THE IMMUNE SYSTEM, 203.
99. Gückel, B., *et al.* (2006). Pre-existing T-cell immunity against mucin-1 in breast cancer patients and healthy volunteers. doi: <https://doi.org/10.1007/s00432-005-0064-6>
100. Zhao, Q., *et al.* (2015). Affinity maturation of T-cell receptor-like antibodies for Wilms tumor 1 peptide greatly enhances therapeutic potential. doi: <https://doi.org/10.1038/leu.2015.125>
101. Oka, Y., *et al.* (2000). Human cytotoxic T-lymphocyte responses specific for peptides of the wild-type Wilms' tumor gene (*WT1*) product.
102. Coelho-dos-Reis, J. G., *et al.* (2016). Co-administration of α -GalCer analog and TLR4 agonist induces robust CD8⁺ T-cell responses to PyCS protein and WT-1 antigen and activates memory-like effector NKT cells. doi: <https://doi.org/10.1016/j.clim.2016.04.014>
103. Lichtenegger, F. S., *et al.* (2018). Targeting LAG-3 and PD-1 to enhance T cell activation by antigen-presenting cells. doi: [10.3389/fimmu.2018.00385](https://doi.org/10.3389/fimmu.2018.00385)
104. Schmied, S., *et al.* (2015). Analysis of the functional WT 1-specific T-cell repertoire in healthy donors reveals a discrepancy between CD 4⁺ and CD 8⁺ memory formation. doi: <https://doi.org/10.1111/imm.12472>
105. Slingluff J., *et al.* (2010). Immunogenicity for CD8⁺ and CD4⁺ T cells of two formulations of an incomplete Freund's adjuvant for multipptide melanoma vaccines. doi: [10.1097/CJI.0b013e3181e311ac](https://doi.org/10.1097/CJI.0b013e3181e311ac)
106. Trefzer, U., *et al.* (2005). Tumour-dendritic hybrid cell vaccination for the treatment of patients with malignant melanoma: immunological effects and clinical results. doi: <https://doi.org/10.1016/j.vaccine.2005.01.081>
107. Kammula, U. S., *et al.* (2008). Use of high throughput qPCR screening to rapidly clone low frequency tumour specific T-cells from peripheral blood for adoptive immunotherapy. doi: <https://doi.org/10.1186/1479-5876-6-60>
108. Border, E. C., *et al.* (2019). Affinity-enhanced T-cell receptors for adoptive T-cell therapy targeting MAGE-A10: strategy for selection of an optimal candidate.
109. Kim, H-D., *et al.* (2021). Implication of CD69⁺ CD103⁺ tissue-resident-like CD8⁺ T cells as a potential immunotherapeutic target for cholangiocarcinoma. doi: <https://doi.org/10.1111/liv.14814>

REFERENCES

110. van der Bruggen, P., *et al.* (1994). A peptide encoded by human gene MAGE-3 and presented by HLA-A2 induces cytolytic T lymphocytes that recognize tumor cells expressing MAGE-3. doi: <https://doi.org/10.1002/eji.1830241218>
111. Choi, C., *et al.* (2005). Enrichment of functional CD8 memory T cells specific for MUC1 in bone marrow of patients with multiple myeloma. doi: <https://doi.org/10.1182/blood-2004-01-0366>
112. Shirakawa, T., *et al.* (2018). Antitumor effect of oral cancer vaccine with Bifidobacterium delivering WT1 protein to gut immune system is superior to WT1 peptide vaccine.
113. Zhou, M., *et al.* (2005). Identification of two naturally presented MAGE antigenic peptides from a patient with hepatocellular carcinoma by mass spectrometry." *Immunology letters* 99.1 (2005): 113-121. doi: <https://doi.org/10.1016/j.imlet.2005.02.007>
114. Kobayashi, T., *et al.* (2003). New MAGE-4 antigenic peptide recognized by cytolytic T lymphocytes on HLA-A1 tumor cells. doi: <https://doi.org/10.1034/j.1399-0039.2003.00123.x>
115. Zhang, Y., *et al.* (2002). A MAGE-A4 peptide presented by HLA-B37 is 392 recognized on human tumors by cytolytic T lymphocytes.
116. Miyahara, Y., *et al.* (2005). Determination of cellularly processed 386 HLA-A2402-restricted novel CTL epitopes derived from 387 two cancer germ line genes, MAGE-A4 and SAGE. doi: [10.1158/1078-0432.CCR-04-2585](https://doi.org/10.1158/1078-0432.CCR-04-2585)
117. Traversari, C., *et al.* (1992). A nonapeptide encoded by human gene MAGE-1 is recognized on HLA-A1 by cytolytic T lymphocytes directed against tumor antigen MZ2-E. doi: <https://doi.org/10.1084/jem.176.5.1453>
118. Gavvovidis, I., *et al.* (2021). Generation of a HuTCR mouse platform-derived MAGE-A1-directed high-affinity TCR with superior potency versus human-derived TCRs.
119. Ottaviani, S., *et al.* (2005). MAGE-1 antigenic peptide recognized by human cytolytic T 246 lymphocytes on HLA-A2 tumor cells. doi: [10.1007/s00262-005-0705-2](https://doi.org/10.1007/s00262-005-0705-2)
120. Wang, Sw., *et al.* (2015). A novel MHC- dextramer assay to identify melanoma antigen-specific CD8+ T cells from solid tumor disaggregates and matched peripheral blood. doi: <https://doi.org/10.1186/2051-1426-3-S2-P109>
121. Goodyear, O., *et al.* (2005). CD8+ T cells specific for cancer germline gene antigens are found in many patients with multiple myeloma, and their frequency correlates with disease burden. doi: [10.1182/blood-2005-02-0563](https://doi.org/10.1182/blood-2005-02-0563)
122. Morgan, R. A., *et al.* (1997). Cancer regression and neurologic toxicity following anti-MAGE-A3 TCR gene therapy. doi: [10.1097/CJI.0b013e3182829903](https://doi.org/10.1097/CJI.0b013e3182829903)
123. Martin, A. D., *et al.* (2021). Re-examination of MAGE-A3 as a T-cell Therapeutic Target. doi: [10.1097/CJI.0000000000000348](https://doi.org/10.1097/CJI.0000000000000348)
124. Pol, J. G., *et al.* (2019). Preclinical evaluation of a MAGE-A3 vaccination utilizing the oncolytic Maraba virus currently in first-in-human trials. doi: <https://doi.org/10.1080/2162402X.2018.1512329>
125. Anderson, L.D., *et al.* (2011). Identification of MAGE-C1 (CT-7) epitopes for T-cell therapy of multiple myeloma. doi: [10.1007/s00262-011-1009-3](https://doi.org/10.1007/s00262-011-1009-3)
126. Tyler, E. M., *et al.* (2013). CT7/MAGE-C1 Expression and Immune Responses Following Allogeneic Transplant for Multiple Myeloma. doi: <https://doi.org/10.1016/j.bbmt.2012.11.116>
127. Casado, J. G., *et al.* (2009). Correlation of effector function with phenotype and cell division after in vitro differentiation of naive MART-1-specific CD8+ T cells. doi: <https://doi.org/10.1093/intimm/dxn123>

REFERENCES

128. Saleh, F. H., et al. (2001). Primary melanoma tumour regression associated with an immune response to the tumour-associated antigen melan-A/MART-1. doi: <https://doi.org/10.1002/ijc.1491>
129. Yamshchikov, G. V., et al. (2001). Evaluation of peptide vaccine immunogenicity in draining lymph nodes and peripheral blood of melanoma patients. doi: [https://doi.org/10.1002/1097-0215\(20010601\)92:5<703::AID-IJC1250>3.0.CO;2-5](https://doi.org/10.1002/1097-0215(20010601)92:5<703::AID-IJC1250>3.0.CO;2-5)
130. Mosse, C. A., et al. (1998). The class I antigen-processing pathway for the membrane protein tyrosinase involves translation in the endoplasmic reticulum and processing in the cytosol. doi: [10.1084/jem.187.1.37](https://doi.org/10.1084/jem.187.1.37)
131. Wölfel, C., et al. (2000). Transporter (TAP)-and proteasome-independent presentation of a melanoma-associated tyrosinase epitope. doi: [https://doi.org/10.1002/1097-0215\(20001101\)88:3<432::AID-IJC16>3.0.CO;2-9](https://doi.org/10.1002/1097-0215(20001101)88:3<432::AID-IJC16>3.0.CO;2-9)
132. Chianese-Bullock, K. A., et al. (2009). Multi-peptide vaccines vialed as peptide mixtures can be stable reagents for use in peptide-based immune therapies. doi: <https://doi.org/10.1016/j.vaccine.2009.01.018>
133. MoreauAubry, A., et al. (1997). Identification of a gp100 epitope recognized by HLA-A3 restricted melanoma infiltrating lymphocytes. <https://doi.org/10.3892/ijo.10.4.841>
134. Wenandy, L., et al. (2008). RhoC a new target for therapeutic vaccination against metastatic cancer. doi: [10.1007/s00262-008-0517-2](https://doi.org/10.1007/s00262-008-0517-2)
135. Reker, S., et al. (2004). Identification of Novel Survivin-Derived CTL Epitopes with Different HLA-A-Restriction Profiles. doi: [10.4161/cbt.3.2.611](https://doi.org/10.4161/cbt.3.2.611)
136. Brewer, K. D., et al. (2018). Unique depot formed by an oil based vaccine facilitates active antigen uptake and provides effective tumour control. doi: <https://doi.org/10.1186/s12929-018-0413-9>
137. Reynolds, S. R., et al. (2000). Identification of HLA-A*03, A*11 and B*07-restricted melanoma-associated peptides that are immunogenic in vivo by vaccine-induced immune response (VIIR) analysis. doi: [10.1016/s0022-1759\(00\)00254-4](https://doi.org/10.1016/s0022-1759(00)00254-4)
138. Nishida, T., et al. (2004). Clinical relevance of a newly identified HLA-A24-restricted minor histocompatibility antigen epitope derived from BCL2A1, ACC-1, in patients receiving HLA genotypically matched unrelated bone marrow transplant. doi: <https://doi.org/10.1111/j.1365-2141.2004.04823.x>
139. Hambach, L., et al. (2005). Immunotherapy of cancer through targeting of minor histocompatibility antigens. doi: <https://doi.org/10.1016/j.coi.2005.01.010>
140. Uemura, H., et al. (2008). Tumor vaccines in renal cell carcinoma. doi: <https://doi.org/10.1007/s00345-008-0246-3>
141. Obara, W., et al. (2012). Cancer peptide vaccine therapy developed from oncoantigens identified through genome-wide expression profile analysis for bladder cancer. doi: <https://doi.org/10.1093/jjco/hys098>
142. Kikuchi, R., et al. (2019). A pilot study of vaccine therapy with multiple glioma oncoantigen/glioma angiogenesis-associated antigen peptides for patients with recurrent/progressive high-grade glioma. doi: <https://doi.org/10.3390/jcm8020263>
143. Daisuke, N., et al. (2013). Peptide vaccines for hepatocellular carcinoma. doi: <https://doi.org/10.4161/hv.22473>
144. Kirkin, A. F., et al. (1998). The immunogenic properties of melanoma-associated antigens recognized by cytotoxic T lymphocytes. doi: <https://doi.org/10.1159/000019050>

REFERENCES

145. Goto, S., *et al.* (2005). Combined immunocell therapy using activated lymphocytes and monocyte-derived dendritic cells for malignant melanoma.
146. Takeuchi, S., *et al.* (2020). Anti-cancer immunotherapy using cancer-derived multiple epitope-peptides cocktail vaccination clinical studies in patients with refractory/persistent disease of uterine cervical cancer and ovarian cancer [phase 2]. doi: <https://doi.org/10.1080/2162402X.2020.1838189>
147. Hasegawa, K., *et al.* (2018). Phase I study of multiple epitope peptide vaccination in patients with recurrent or persistent cervical cancer. doi: <https://doi.org/10.1097/CJI.0000000000000214>
148. Osawa, R., *et al.* (2012). Identification of HLA-A24-restricted novel T Cell epitope peptides derived from P-cadherin and kinesin family member 20A. doi: <https://doi.org/10.1155/2012/848042>
149. Suzuki, H., *et al.* (2013). Multiple therapeutic peptide vaccines consisting of combined novel cancer testis antigens and anti-angiogenic peptides for patients with non-small cell lung cancer. doi: <https://doi.org/10.1186/1479-5876-11-97>
153. Jiang, X., *et al.* (2007). Induction of cytotoxic T lymphocytes specific to malignant glioma using T2 cells pulsed with HLA-A2-restricted interleukin-13 receptor $\alpha 2$ peptide in vitro. doi: <https://doi.org/10.1111/j.1745-7270.2007.00331.x>
154. Termine, J. D. (1988). Non-collagen proteins in bone." Cell and molecular biology of vertebrate hard tissues.
155. Yoshitake, Y., *et al.* (2015). Phase II clinical trial of multiple peptide vaccination for advanced head and neck cancer patients revealed induction of immune responses and improved OS. doi: [10.1158/1078-0432.CCR-14-0202](https://doi.org/10.1158/1078-0432.CCR-14-0202)
156. Suda, T., *et al.* (2007). Identification of human leukocyte antigen-A24-restricted epitope peptides derived from gene products upregulated in lung and esophageal cancers as novel targets for immunotherapy. doi: <https://doi.org/10.1111/j.1349-7006.2007.00603.x>
157. Imamura, Y., *et al.* (2016). Generation of Large Numbers of Antigen-Expressing Human Dendritic Cells Using CD14-ML Technology. doi: <https://doi.org/10.1371/journal.pone.0152384>
158. Suzuki, H., *et al.* (2013). Multiple therapeutic peptide vaccines consisting of combined novel cancer testis antigens and anti-angiogenic peptides for patients with non-small cell lung cancer. doi: <https://doi.org/10.1186/1479-5876-11-97>
159. Takeda, K., *et al.* (2018). Quantitative T-cell repertoire analysis of peripheral blood mononuclear cells from lung cancer patients following long-term cancer peptide vaccination. doi: <https://doi.org/10.1007/s00262-018-2152-x>
160. Okuyama, R., *et al.* (2013). Immunological responses to a multi-peptide vaccine targeting cancer-testis antigens and VEGFRs in advanced pancreatic cancer patients. doi: <https://doi.org/10.4161/onci.27010>
161. Aruga, A., *et al.* (2014). Phase I clinical trial of multiple-peptide vaccination for patients with advanced biliary tract cancer. doi: <https://doi.org/10.1186/1479-5876-12-61>
162. Wang, Y., *et al.* (2006). Prostate-specific antigen-reactive cytotoxic T lymphocyte precursors in colon cancer patients. doi: <https://doi.org/10.3892/or.15.2.317>
163. Harada, M., *et al.* (2003). Prostate-specific antigen-derived epitopes capable of inducing cellular and humoral responses in HLA-A24+ prostate cancer patients. doi: <https://doi.org/10.1002/pros.10280>

REFERENCES

164. Hirohashi, Y., *et al.* (2002). An HLA-A24-restricted cytotoxic T lymphocyte epitope of a tumor-associated protein, survivin. *Clin Cancer Res.*
165. Tsuruma, T., *et al.* (2004). Phase I clinical study of anti-apoptosis protein, survivin-derived peptide vaccine therapy for patients with advanced or recurrent colorectal cancer. doi: <https://doi.org/10.1186/1479-5876-2-19>
166. Miyazaki, A., H., *et al.* (2011). Phase I clinical trial of survivin-derived peptide vaccine therapy for patients with advanced or recurrent oral cancer. doi: <https://doi.org/10.1111/j.1349-7006.2010.01789.x>
167. Suda, T., *et al.* (2007). Identification of human leukocyte antigen-A24-restricted epitope peptides derived from gene products upregulated in lung and esophageal cancers as novel targets for immunotherapy. doi: <https://doi.org/10.1111/j.1349-7006.2007.00603.x>
168. Aruga, A., *et al.* (2013). Vaccination of biliary tract cancer patients with four peptides derived from cancer-testis antigens. doi: [10.4161/onci.24882](https://doi.org/10.4161/onci.24882)
169. Tsuboi, A., *et al.* (2002). Enhanced induction of human WT1-specific cytotoxic T lymphocytes with a 9-mer WT1 peptide modified at HLA-A*2402-binding residues. doi: <https://doi.org/10.1007/s00262-002-0328-9>
170. Sawada, A., *et al.* (2016). Feasibility of Cancer Immunotherapy with WT1 Peptide Vaccination for Solid and Hematological Malignancies in Children. doi: <https://doi.org/10.1002/pbc.25792>
171. Matsuda T, *et al.* (2018). Pilot study of WT1 peptide-pulsed dendritic cell vaccination with docetaxel in esophageal cancer.
172. Kang, X., *et al.* (1995). Identification of a tyrosinase epitope recognized by HLA-A24-restricted, tumor-infiltrating lymphocytes. *The J. of Immunol.*
173. Morck, T. A., *et al.* (2003). Mesothelin-specific CD8+ T cell responses provide evidence of in vivo cross-priming by antigen-presenting cells in vaccinated pancreatic cancer patients. doi: <https://doi.org/10.1084/jem.20031435>
174. Matsuzaki, J., *et al.* (2008). Recognition of naturally processed and ovarian cancer reactive CD8+ T cell epitopes within a promiscuous HLA class II T-helper region of NY-ESO-1. doi: <https://doi.org/10.1007/s00262-008-0450-4>
175. Stager, E. H., *et al.* (2004). Induction of CAMEL/NY-ESO-ORF2-specific CD8+ T cells upon stimulation with dendritic cells infected with a modified Ad5 vector expressing a chimeric Ad5/35 fiber. doi: [10.1038/sj.cgt.7700674](https://doi.org/10.1038/sj.cgt.7700674)
176. Wieger, J., *et al.* (2009). Myeloid leukemic progenitor cells can be specifically targeted by minor histocompatibility antigen LRH-1-reactive cytotoxic T cells. doi: <https://doi.org/10.1182/blood-2008-04-153825>
177. Hobo, W., *et al.* (2013). Improving dendritic cell vaccine immunogenicity by silencing PD-1 ligands using siRNA-lipid nanoparticles combined with antigen mRNA electroporation. doi: <https://doi.org/10.1007/s00262-012-1334-1>
178. Luiten, R., and Van Der Bruggen, P., *et al.* (2000). A MAGE-A1 peptide is recognized on HLA-B7 human tumors by cytolytic T lymphocytes. doi: <https://doi.org/10.1034/j.1399-0039.2000.550206.x>
179. Borobova, E. A., *et al.* (2018). Design of Artificial Immunogens Containing Melanoma-associated T-cell Epitopes. doi: <https://doi.org/10.2174/1566523218666181113112829>

REFERENCES

180. Benlalam, H., *et al.* (2003). Identification of five new HLA-B*3501-restricted epitopes derived from common melanoma-associated antigens, spontaneously recognized by tumor-infiltrating lymphocytes. doi: [10.4049/jimmunol.171.11.6283](https://doi.org/10.4049/jimmunol.171.11.6283)
181. Postow, M. A., *et al.* (2014). Markers for anti-cytotoxic T-lymphocyte antigen 4 (CTLA-4) therapy in melanoma. doi: https://doi.org/10.1007/978-1-62703-727-3_6
182. Lin, Y., *et al.* (2009). Optimization and validation of a robust human T-cell culture method for monitoring phenotypic and polyfunctional antigen-specific CD4 and CD8 T-cell responses. doi: [10.3109/14653240903136987](https://doi.org/10.3109/14653240903136987)
183. Engelstein, R., *et al.* (2016). Immune monitoring of patients treated with a whole-cell melanoma vaccine engineered to express 4-1BBL. doi: [10.1097/cji.0000000000000138](https://doi.org/10.1097/cji.0000000000000138)
184. Kitano, S., *et al.* (2013). Enhancement of Tumor-Reactive Cytotoxic CD4. doi: [10.1158/2326-6066.CIR-13-0068](https://doi.org/10.1158/2326-6066.CIR-13-0068)
185. Jäger, E., *et al.* (2002) Identification of a naturally processed NY-ESO-1 peptide recognized by CD8+ T cells in the context of HLA-B51. PMID: 12747757
186. Zeng G, *et al.* (2002). Generation of NY-ESO-1-specific CD4+ and CD8+ T cells by a single peptide with dual MHC class I and class II specificities: a new strategy for vaccine design. PMID: [PMC2562286](https://pubmed.ncbi.nlm.nih.gov/12562286/)
187. Hunder, N. N., *et al.* (2008). Treatment of metastatic melanoma with autologous CD4+ T cells against NY-ESO-1.
188. Birkholz, K., *et al.* (2010). Targeting of DEC-205 on human dendritic cells results in efficient MHC class II-restricted antigen presentation. doi: <https://doi.org/10.1182/blood-2010-02-268425>
189. Schultz, E. S., *et al.* (2000). A MAGE-A3 peptide presented by HLA-DP4 is recognized on tumor cells by CD4+ cytolytic T lymphocytes. Cancer Res.
190. Yao, X., *et al.* (2016). Isolation and Characterization of an HLA-DPB1*04: 01-restricted MAGE-A3 T-Cell Receptor for Cancer Immunotherapy. doi: [10.1097/CJI.0000000000000123](https://doi.org/10.1097/CJI.0000000000000123)
191. Lowe, D. B., *et al.* (2013). Combination Therapies to Improve Delivery of Protective T Cells into the Melanoma Microenvironment. Melanoma - From Early Detection to Treatment. doi: 10.5772/53634. Available from: <https://www.intechopen.com/chapters/42292>
192. Houghton C. S. B., *et al.* (2007). Immunological validation of the EpiOptimizer program for streamlined design of heteroclitic epitopes. doi: <https://doi.org/10.1016/j.vaccine.2007.05.008>