

**Reagent provided**

TCR monomer is uniquely identified by its account product number, e.g., MTCRU0026. Soluble TCR Monomer consists of a TCR alpha and beta chain. The TCR monomer is mono-biotinylated.

The product is provided in buffer containing 20 mM Tris-HCl, 50 mM NaCl, pH 8.0.

This is a custom-made product based on TCR alpha and beta chain sequences provided by customer.

For research use only. Not for use in diagnostic or therapeutic procedures

**Sizes**

See amount and concentration on label

**Storage**

Store at -80°C. Avoid repeat freeze-thaw cycles.

**Precautions**

As with any product derived from biological sources, proper handling procedures should be used.

For professional users.

**Symbols**

See [www.immudex.com/symbols](http://www.immudex.com/symbols)

**Technical support**

E-mail: [customer@immudex.com](mailto:customer@immudex.com)

Telephone: +45 3110 9292 (Denmark)

**Manufacturer**

Immudex, Bredevej 2A, DK-2830 Virum, Denmark